



Aug 12,
2019

untitled protocol V.1

Khairil Asnan¹, Mohd. Sadzan Asram bin Aliming², Aprialdy Idrus², Diana Eka Pratiwi²

¹Biology Department, Universitas Negeri Makassar, ²Chemistry Department, Universitas Negeri Makassar

In Development

[dx.doi.org/10.17504/protocols.io.6fyhbpw](https://doi.org/10.17504/protocols.io.6fyhbpw)



Khairil Asnan



ABSTRACT

Cattle's blood component is compressed using centrifuge (5000 rpm for 10 minutes), then supernatant (erythrocytes) collected. Erythrocytes with a ratio of 3:1 (serum:supernatant), treated with initial mixing each with HCl (0.5 M) NaOH (0.5 M), and 30% H₂O₂ with an initial ratio of 1:1 (v/v). Each variable then treated by mixing SnCl₂ (1 M, 1:1 v/v), oleic acid (2:1 v/v), and O_{2(g)} (0.5 L/m, 60s). The best sample then purified by heat about 800°C for 120 minutes. The highest calculated total iron mass is 240000 µg/100 ml (NaOH + oleic acid). Sample purification has increased the iron concentration up to 46.30% (m/m%).

GUIDELINES

Beware to react the blood with peroxide because it will produces bubbles

MATERIALS

NAME	CATALOG #	VENDOR
Sodium Hydroxide	View	Sigma Aldrich

STEPS MATERIALS

NAME	CATALOG #	VENDOR
Sodium Hydroxide	View	Sigma Aldrich

MATERIALS TEXT

[M] 0.5 Molarity (M)

37 °C

00:30:00



Sodium Hydroxide




by Sigma Aldrich

[View](#)

10 rpm Mix with sample for 20 seconds

Sample preparation

1

Compress the whole blood sample from cattle.  6 ml  37 °C  00:10:00








Mini-centrifuge
Centrifuge

Fisher S67601B [↗](#)

Any standard mini centrifuge with adapters for different tube sizes will suffice



 5000 rpm

2 Mix  2 ml blood sample sample with  2 ml NaOH  37 °C  00:00:30
 0.5 Molarity (M)







Sodium Hydroxide
by Sigma Aldrich
[View](#)



dark green solution with strong odor

 10 rpm after we rest the mixed solution for 30 seconds

- 3 Chelate reaction by adding  **4 ml Oleic acid** into the previous mixed sample solution then let it sit for  **00:00:30** . After that, mixed them by  **10 rpm** for  **00:00:30** .



Centrifuge

Benchtop Centrifuge

Eppendorf 5405000441 

Any benchtop centrifuge will suffice



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited