



Feb 01, 2020

## Initial Rapid Pathology Assessment of Kidney Tissue

Jamie Allen<sup>1</sup>, Elizabeth Neumann<sup>1</sup>, Maya Brewer<sup>1</sup>, Haichun Yang<sup>1</sup>, Jeff Spraggins<sup>1</sup>, Danielle Gutierrez<sup>1</sup>

<sup>1</sup>Vanderbilt University

1 Works for me dx.doi.org/10.17504/protocols.io.9dph25n

VU Biomolecular Multimodal Imaging Center

Human BioMolecular Atlas Program (HuBMAP) Method Development Community



### ABSTRACT

#### Scope:

Prepare formalin fixed tissue with freezing fresh tissue protocol (dx.doi.org/10.17504/protocols.io.6wghfbw).

This protocol provides the steps to prepare formalin fixed human kidney tissue and perform histology assessment for normalcy.

#### Expected Outcome:

Pathological assessment of kidney tissue for use in LC-MS/MS and imaging.

- 1 Once fixed (dx.doi.org/10.17504/protocols.io.6wghfbw) take tissue to pathology Core for paraffin emedding.
- 2 Section samples at 5 µm.
- 3 PAS stain tissue sections (dx.doi.org/10.17504/protocols.io.4qngvve).
- 4 Scan slides with brightfield scanner (Leica) and save as .tiff or .jpg
- 5 Place saved images on QuPath for analysis.
- 6 Assess and record the following information for each tissue:  
(%): cortex and medulla  
(Yes/No): Pyramid Presence, Autolysis, Non-Renal disease (*i.e.* cancer)  
(0=None – 3=Severe): Glomerular disease, Tubulointerstitial disease
- 7 Based on FFPE assessment, kidneys will be used accordingly:  
  
LC-MS/MS normalcy:  
low or no glomerular disease, and 75% or more cortex for normal  
100% tumor for diseased  
  
3D Imaging:  
50:50 cortex:medulla  
low or no glomerular disease  
no renal disease



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited