



## Human RBM3 ELISA

[Lisa-Maria Rosenthal](#)<sup>1</sup>, Giang Tong<sup>1</sup>, Katharina Schmitt<sup>1</sup>

<sup>1</sup>Deutsches Herzzentrum Berlin

[dx.doi.org/10.17504/protocols.io.3p5gmq6](https://doi.org/10.17504/protocols.io.3p5gmq6)

Working

 [Lisa-Maria Rosenthal](#)  
Deutsches Herzzentrum Berlin 

## MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Human RBM3 ELISA Kit

not available (prepared for cust

- 1 Prepare RBM3 standards in RBM3 Dilution Buffer: 0, 31.25, 62.5, 125, 250, 500, 1000 and 2000 pg/ml.
- 2 Dilute the serum samples 1:10 with RBM3-SAMPLE-BUF.
- 3 Pipette 50µl/well RBM3-Standards (0-2000 pg/ml) and 50µl/well diluted samples in duplicates onto 96-well microplate.
- 4 Pipette 50µl/well detection conjugate – diluted RBM3-HRP conjugate onto 96-well microplate.
- 5 Shake for 5 min and incubate for 120 min at 20-25 °C.
- 6 Wash 3 times with 1:10 diluted Wash solution.
- 7 Pipette 50µl/well Tetramethylbenzidin solution and incubate for 30 min at 20-25 °C in the dark on a shaker.
- 8 Pipette 50µl/well Stop Reagent.
- 9 Measure absorbance at 450 nm using a plate reader.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited