

Biolistic Transformation of *Pseudo-nitzschia multistriata* and *Pseudo-nitzschia arenysensis*

Valeria Sabatino, Monia Teresa Russo, Shrikant Patil, Giuliana d'Ippolito, Angelo Fontana, Maria Immacolata Ferrante

Abstract

This protocol allows transformation of two species of *Pseudo-nitzschia* through biolistic method.

With this method it is possible to insert one (transformation) or two (co-transformation) plasmids in the cells.

Citation: Valeria Sabatino, Monia Teresa Russo, Shrikant Patil, Giuliana d'Ippolito, Angelo Fontana, Maria Immacolata Ferrante Biolistic Transformation of *Pseudo-nitzschia multistriata* and *Pseudo-nitzschia arenysensis*. **protocols.io** dx.doi.org/10.17504/protocols.io.hp5b5q6

Published: 19 Apr 2017

Protocol

To be completed in advance

Step 1.

Prepare plasmid with your gene of interest at the concentration of $1 \mu\text{g ml}^{-1}$.

To be completed in advance

Step 2.

Grow wild-type cells in f/2 medium up to mid log-phase.

To be completed in advance

Step 3.

Prepare f/2 plates (0.4% agarose in f/2 medium) and leave to dry under hood.

The day before bombarding

Step 4.

Count the cells. For each plate between $3\text{-}5 \times 10^7$ cells are needed.

The day before bombarding

Step 5.

Centrifuge the cells at 2465 g and resuspend the pellet in a very small volume of f/2 medium.

The day before bombarding

Step 6.

Spread the cells onto agarose plates. Leave to dry under hood.

The day of the transformation

Step 7.

DNA precipitation is carried out using Tungsten M-10 microcarrier 0.7 μm particles (BioRad).

The day of the transformation

Step 8.

Dilute the tungsten particles in a final concentration of 3 mg in 50 μl of water for plate.

The day of the transformation

Step 9.

Add to particles 5 μl of plasmid (1 $\mu\text{g } \mu\text{l}^{-1}$) (for the transformation) or 3 μl of each plasmid (1 $\mu\text{g } \mu\text{l}^{-1}$) (for the co-transformation), 50 μl 2.5 M CaCl_2 , and 20 μl of 0.1 M spermidine. Continue vortexing for 3 min, then microfuge for 10 sec and remove supernatant. Wash with 250 μl 100% ethanol, microfuge again, then resuspend in 60 μl ethanol. Keep in ice.

The day of the transformation

Step 10.

Wash rupture disks and macrocarriers in ethanol and allow to dry.

Bombardment

Step 11.

Bombardment is carried out in a Biorad PDS-1000/He microprojectile accelerator. Shooting is carried out using 1,550 PSI rupture disks. Plates for transformation are placed in the middle slot immediately under the macrocarrier.

Bombardment

Step 12.

10 μl aliquot of the prepared particles were transferred to the centre of each macrocarrier. Allow to dry.

Bombardment

Step 13.

Keep the plates in incubator overnight.

The day after the bombardment

Step 14.

Resuspend the cells (each plate) in 200 ml of f/2 medium with antibiotic. *Pseudo-nitzschia* cells grow better in liquid medium.

Step 15.

Usually untransformed *Pseudo-nitzschia* cells are dead at 10-12 days post-application of antibiotic.