



Oct 24, 2019

Miracle Prep for Plasmid Isolation

Blake Flood¹¹University of Chicago

1

Works for me

dx.doi.org/10.17504/protocols.io.8p2hvhqe



Blake Flood ⚡

ABSTRACT

Based on:

<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0160509>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Pronobis MI, Deutch N, Peifer M (2016) The Miraprep: A Protocol that Uses a Miniprep Kit and Provides Maxiprep Yields. PLoS ONE 11(8): e0160509. <https://doi.org/10.1371/journal.pone.0160509>

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Sterile deionized H2O		
LB Broth		
Buffer P1	19051	Qiagen
Buffer P2	19052	Qiagen
Buffer N3	19064	Qiagen
96% EtOH		
Econospin Mini Spin Column 250ct	1920-250	Epoch Life Science
Buffer PB	19066	Qiagen
Buffer PE	19065	Qiagen

STEPS MATERIALS




NAME ▾	CATALOG # ▾	VENDOR ▾
Buffer P1	19051	Qiagen
Buffer P2	19052	Qiagen
Buffer N3	19064	Qiagen
Buffer PB	19066	Qiagen

Grow bacteria containing plasmid

- 1 Pick colony from plate. Inoculate 50 mL LB broth per vector.

2 Allow to grow overnight at  **37 °C** .

Isolate Plasmid

3 Spin culture  **00:10:00**  **4000 x g**  **4 °C**

4 Remove supt. Resuspend cells in  **2 ml** Buffer P1.

5  **2 ml** Buffer P2 (invert 3-4 times)



Buffer P1

by Qiagen

Catalog #: [19051](#)

6 Incubate  **00:03:00** @  **Room temperature**



Buffer P2

by Qiagen

Catalog #: [19052](#)

7  **2 ml** Buffer N3 (invert 3-4 times)



Buffer N3


by Qiagen


Catalog #: [19064](#)

8 Distribute lysate into  **4 x** 1.5mL tubes

9 Centrifuge  **00:10:00** @  **13200 x g**

10 Collect supt, add 1 volume 96% EtOH, mix

11 Load on  **4 x** "QIAprep 2.0 spin columns" or any generic silica membrane mini spin column (Econospin).

12 Spin  **00:01:00** @  **17900 x g**

13 Re-load and re-spin until everything is loaded.

14  **500 µl** Buffer PB to each column



Buffer PB

by QIAGEN

Catalog #: 19066

15 Spin ⌚ 00:01:00 @ 🌀 17900 x g

16 750 uL Buffer PE

17 Spin ⌚ 00:01:00 @ 🌀 17900 x g

18 Spin again to dry ⌚ 00:01:00 @ 🌀 17900 x g

19 Elute with 30-50 uL H₂O

20 (Optional) Re-load H₂O on columns and re-spin ⌚ 00:01:00 @ 🌀 17900 x g to increase yield.

step case

If need higher yield

step case

If need higher yield



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited