



Apr 16, 2019

Working

Cell growth assay [↗](#)

PLOS One

Kiichi Hirota¹, Yoshiyuki Matsuo¹¹Kansai Medical University[dx.doi.org/10.17504/protocols.io.v64e9gw](https://doi.org/10.17504/protocols.io.v64e9gw)

Yoshiyuki Matsuo

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0215072>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Sumi C, Matsuo Y, Kusunoki M, Shoji T, Uba T, Iwai T, Bono H, Hirota K (2019) Cancerous phenotypes associated with hypoxia-inducible factors are not influenced by the volatile anesthetic isoflurane in renal cell carcinoma. PLoS ONE 14(4): e0215072. doi: [10.1371/journal.pone.0215072](https://doi.org/10.1371/journal.pone.0215072)

BEFORE STARTING

Cell growth was assessed using a CellTiter 96® Aqueous One Solution Cell Proliferation Assay (Promega, Madison, WI, USA). This assay measured the reduction of the tetrazolium compound MTS (3-[4,5-dimethyl-2-yl]-5-[3-carboxymethoxyphenyl]-2-[4-sulphophenyl]-2H-tetrazolium, inner salt).

- 1 Cells were seeded into 96-well plates and cultivated for indicated time periods
- 2 20 µl of CellTiter 96® Aqueous One Solution Reagent was added to each well.
- 3 The plates were incubated at 37 °C for 1 h prior to measuring the absorbance of each sample using an iMark™ Microplate Reader (BIO-RAD, Hercules, CA, USA) at a wavelength of 490 nm.
- 4 Cell viability was then calculated by comparing the absorbance of the treated cells with that of the control cells, with the latter defined as 100%.

All samples were assayed in triplicate or quadruplicate for each experiment.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited