

NEBNext Ultra II FS DNA Module E7810

New England Biolabs¹

¹New England Biolabs

dx.doi.org/10.17504/protocols.io.n34dgqw

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com

ABSTRACT

The NEBNext Ultra II FS DNA Module contains the enzymes and buffers required to convert a broad range of input amounts of intact DNA into fragmented DNA with 5' phosphorylated 3' dA-tailed ends. The module is optimized for use with the NEBNext Ultra II Ligation Module (NEB #E7595) and with the NEBNext Ultra II Q5 Master Mix (NEB #M0544) if amplification is required. The fast, user-friendly workflow has minimal hands on time.

Note: The Ultra II FS Module is not compatible with bisulfite conversion workflows



manuale7810 1.0
9:17.pdf

PROTOCOL STATUS

Working

GUIDELINES

The NEBNext Ultra II FS DNA Module is Designed for use with the Following:

NEBNext Ultra II Ligation Module (NEB #E7595)

NEBNext Ultra II Q5® Master Mix (NEB #M0544)

NEBNext Singleplex or Multiplex Oligos for Illumina®




(NEB #E7350, #E7335, #E7500, #E6609, #E7710, #E7730 or #E7600)

MATERIALS

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|-------------------------------------|-------------|---------------------|
| TE Buffer (1X) | E7808 | New England Biolabs |
| NEBNext Ultra II FS Reaction Buffer | E7807 | New England Biolabs |
| NEBNext Ultra II FS Enzyme Mix | E7806 | New England Biolabs |
| Vortex | View | |
| Microcentrifuge | View | |
| 0.2 ml thin wall PCR tubes | View | |
| PCR Machine | View | |

STEPS MATERIALS

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|-------------------------------------|-------------|---------------------|
| NEBNext Ultra II FS Reaction Buffer | E7807 | New England Biolabs |

| NAME ▾ | CATALOG # ▾ | VENDOR ▾ |
|-------------------------------------|-------------|---|
| NEBNext Ultra II FS Enzyme Mix | E7806 |  New England Biolabs |
| NEBNext Ultra II FS Reaction Buffer | E7807 |  New England Biolabs |
| NEBNext Ultra II FS Enzyme Mix | E7806 |  New England Biolabs |

BEFORE STARTING


Starting Material: 100 pg–500 ng purified, genomic DNA. We recommend that the DNA be in 1X TE (10 mM Tris pH 8.0, 1 mM EDTA), however, 10 mM Tris pH 7.5–8, low EDTA TE or H₂O are also acceptable. If the input DNA is less than 26 µl, add TE (provided) to a final volume of 26 µl.

Fragmentation/End Prep

- 1 Fragmentation occurs during the 37°C incubation step. Use the chart below to determine the incubation time required to generate the desired fragment sizes. Incubation time may need to be optimized for individual samples. See Figure 1 for a typical fragmentation pattern.

| Fragmentation Size | Incubation @ 37°C | Optimization |
|--------------------|-------------------|--------------|
| 100 bp–250 bp | 30 min | 30–40 min |
| 150 bp–350 bp | 20 min | 20–30 min |
| 200 bp–450 bp | 15 min | 15–20 min |
| 300 bp–700 bp | 10 min | 5–15 min |
| 500 bp–1 kb | 5 min | 5–10 min |

Figure 1: Example of size distribution on a Bioanalyzer®. Human DNA (NA19240) was fragmented for 5–40 min.



- 2 Ensure that the Ultra II FS Reaction Buffer is completely thawed. If a precipitate is seen in the buffer, pipette up and down several times to break it up, and quickly vortex to mix. Place on ice until use.



NEBNext Ultra II FS Reaction Buffer
by [New England Biolabs](#)
Catalog #: [E7807](#)

- 3 Vortex the Ultra II FS Enzyme Mix 5-8 seconds prior to use and place on ice.

NOTE

It is important to vortex the enzyme mix prior to use for optimal performance.



NEBNext Ultra II FS Enzyme Mix

by [New England Biolabs](#)

Catalog #: [E7806](#)

- 4 Add the following components to a 0.2 ml thin wall PCR tube on ice:

| Component | Volume per One Library |
|-------------------------------------|------------------------|
| DNA | 26 μ l |
| NEBNext Ultra II FS Reaction Buffer | 7 μ l |
| NEBNext Ultra II FS Enzyme Mix | 2 μ l |
| Total Volume | 35 μ l |



NEBNext Ultra II FS Reaction Buffer

by [New England Biolabs](#)

Catalog #: [E7807](#)



NEBNext Ultra II FS Enzyme Mix

by [New England Biolabs](#)

Catalog #: [E7806](#)

- 5 Vortex the reaction for  00:00:05 and briefly spin in a microcentrifuge.

- 6 In a thermocycler, with the heated lid set to 75°C, run the following program:

5–30 min @ 37°C

30 min @ 65°C

Hold @ 4°C

NOTE

If necessary, samples can be stored at –20°C; however, a slight loss in yield (~20%) may be observed. We recommend continuing with adaptor ligation using the NEBNext Ultra II Ligation Module (NEB #E7595) before stopping.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited