

Field sampling of root-associated microbes for DNA/RNA extraction

Roey Angel

Abstract

This protocol describes a procedure for sampling plant roots in the field for future DNA and RNA extraction for microbiome analysis. The protocol is deliberately designed to be simple and requires no electronic equipment. Root samples are preserved in LifeGuard Soil Preservation Solution for protecting against nucleic acid degradation.

Citation: Roey Angel Field sampling of root-associated microbes for DNA/RNA extraction. **protocols.io**

dx.doi.org/10.17504/protocols.io.pjkdckw

Published: 18 Apr 2018

Before start

Clean spatulas using 70% ethanol

Materials

Micro-spatula set [AT16.1](#) by [Carl Roth](#)

LifeGuard Soil Preservation Solution [12868-100](#) by [Qiagen](#)

Scissors [HCT7.1](#) by [Carl Roth](#)

Technical-grade ethanol (70%) [T913.1](#) by [Carl Roth](#)

Paper towels [Y03.1](#) by [Carl Roth](#)

Microcentrifuge tubes 2 ml [CK06.1](#) by [Carl Roth](#)

Garden trowel [View](#) by [Amazon](#)

Disposable pasteur pipettes [EA61.1](#) by [Carl Roth](#)

Tweezers set [PX40.1](#) by [Carl Roth](#)

Cooling box [AA46.1](#) by [Carl Roth](#)

Cooling packs [E447.1](#) by [Carl Roth](#)

Protocol

Step 1.

Sample a triplicate of plant individuals spaced out a few metres apart from each other

Make sure you are really sampling individual plants and not offshoots of the same plant

Step 2.

Using a garden trowel, carefully dig out the plants while keeping the root system intact (as much as possible, of course)



REAGENTS

Garden trowel [View](#) by [Amazon](#)



Step 3.

While holding the plant by the shoot, shake the root system hard enough so that all loose soil is removed from it.

Take care to damage the plant as little as possible

You can use a spatula to remove large soil aggregates that are attached to the roots

G
M
T

☐ Detect language  English 

Text-to-speech function is limited to 200 characters



[Options](#) : [History](#) : [Feedback](#) : [Donate](#)

Close



REAGENTS

Micro-spatula set [AT16.1](#) by [Carl Roth](#)

Step 4.

From the remaining root system (plus soil particles plus attached to the roots), trim a 'representative sample' of roots using scissors or scalpel

It is usually best to trim the roots onto a piece of paper towel



REAGENTS

Scissors [HCT7.1](#) by [Carl Roth](#)

Paper towels [Y03.1](#) by [Carl Roth](#)

Step 5.

Cut the trimmed out roots a little so that they fit into a 2.0 ml tube



REAGENTS

Scissors [HCT7.1](#) by [Carl Roth](#)

Paper towels [Y03.1](#) by [Carl Roth](#)

Step 6.

Place about 2-3 g of that cut out sample into a 2.0 ml tube



AMOUNT

2 g Additional info:



REAGENTS

Tweezers set [PX40.1](#) by [Carl Roth](#)

Microcentrifuge tubes 2 ml [CK06.1](#) by [Carl Roth](#)



NOTES

Roey Angel 18 Apr 2018

The root tissue should make up at least half or more of the mass, while the remaining attached soil should make up the rest

Step 7.

Press the sample a little into the bottom of the tube to decrease its volume



REAGENTS

Micro-spatula set [AT16.1](#) by [Carl Roth](#)



NOTES

Roey Angel 18 Apr 2018

Make sure the roots do not take up more than 1/2-2/3 of the volume

It is of course possible to split each sample into several separate tubes, depending on the specific type of roots, and submerge each with LifeGuard solution

Step 8.

Add as much LifeGuard solution so that the sample is submerged in about twice of its volume (about 1.0 - 1.5 ml)

Best is to use a disposable Pasteur pipette for dispensing the solution



AMOUNT

1.5 ml Additional info:



REAGENTS

LifeGuard Soil Preservation Solution [12868-100](#) by [Qiagen](#)

Disposable pasteur pipettes [EA61.1](#) by [Carl Roth](#)

Step 9.

Place the tubes in cooling (around 4 °C) and keep them cooled until you reach the lab. The solution will protect nucleic acids even at room temperature for several days, but cooling is preferred

TEMPERATURE

4 °C Additional info:



REAGENTS

Cooling box [AA46.1](#) by [Carl Roth](#)

Cooling packs [E447.1](#) by [Carl Roth](#)

Step 10.

In the lab, store the samples in a freezer (-20 – -80 °C)

TEMPERATURE

-20 °C Additional info: or -80 °C