

# Genomic DNA extraction and PCR

Mark Dewitt & Julia Wong

## Abstract


**Citation:** Mark Dewitt & Julia Wong Genomic DNA extraction and PCR. **protocols.io**

[dx.doi.org/10.17504/protocols.io.dkt4wm](https://dx.doi.org/10.17504/protocols.io.dkt4wm)


**Published:** 12 Aug 2015

## Guidelines

Example of two PCR'd gDNA extracts, one from a nucleofected cell and one from an unzapped control, saw heteroduplexing in the edited sample and tried to get rid of it by reconditioning and/or adding more polymerase.



2% agarose, HBB region PCR'd from K562 gDNA



RC = 10-fold dilution in master mix (same composition), additional 3-cycle PCR  
 E = add 2x amount of enzyme to original master mix  
 RC + E = 10-fold dilution in master mix with double amount of enzyme, additional 3-cycle PCR

## Protocol

### Extraction

#### Step 1.

Count cells

### Extraction

#### Step 2.

Determine volume of cells and volume of QuickExtract DNA extraction solution needed to achieve concentration of 500K cells per 200  $\mu$ L QE solution.

### Extraction

#### Step 3.

Spin down appropriate volume of cells: 300xg, 10 minutes

⌚ DURATION

00:10:00

### Extraction

#### Step 4.

Resuspend in appropriate volume of QE solution, transfer to 1.5mL tubes

### Extraction

#### Step 5.

Pre-heat heating block to 65°C

### Extraction

#### Step 6.

Vortex samples thoroughly for 15 seconds each.

⌚ DURATION

00:00:15

### Extraction

#### Step 7.

Heat samples on block (65°C) for 6 minutes

## DURATION

00:06:00

### Extraction

#### Step 8.

Take samples off block, pre-heat block to 98°C (It's very important for block to reach this temperature! Not enough heat can lead to a failed extraction.)

### Extraction

#### Step 9.

Vortex samples again, 15 seconds each.

## DURATION

00:00:15

### Extraction

#### Step 10.

Heat samples on block (98°C) for 2 minutes.

## DURATION

00:02:00

### Initial PCR

#### Step 11.

Mix for 1rxn

uL	component
29	H <sub>2</sub> O
10	5x Q5 Reaction Buffer (NEB)
1	10mM dNTPs
2.5	10uM Forward Primer
2.5	10uM Reverse Primer
0.5	Q5 DNA polymerase (NEB 2U/uL)
5	extract DNA

run PCR:

1.98C 30s

2.98C 10s

3.66C 30s

4. 72C 20s

5. #2-4, 27 cycles

6.72C 2min\*

\*This final extension time may not be necessary for our short amplicons, and might be contributing to unwanted re-annealing and heteroduplexing.

The products of this initial PCR often contain heteroduplexes that can interfere with downstream analysis. To get rid of heteroduplexes, we use a "reconditioning" step (see next.)

### Reconditioning PCR

### Step 12.

The products of this initial PCR often contain heteroduplexes that can interfere with downstream analysis. To get rid of heteroduplexes, we use a "reconditioning" step (see next.)

### Reconditioning PCR

### Step 13.

Reconditioning PCR

This is essentially a 10-fold dilution of the initial PCR products in more master mix.

Mix for 1rxn

uL	component
29	H2O
10	5x Q5 Reaction Buffer (NEB)
1	10mM dNTPs
2.5	10uM Forward Primer
2.5	10uM Reverse Primer
1.0*	Q5 DNA polymerase (NEB 2U/uL)
5	INITIAL PCR PRODUCT

\*0.5uL is fine too, but we have seen that using double the concentration of polymerase in the reconditioning master mix helps to further eliminate heteroduplexing.

run PCR:

1.98C 30s

2.98C 10s

3.66C 30s

4. 72C 20s

5. #2-4, 27 cycles

no final extension

### Step 14.

Gel: Add 1-2  $\mu$ L 6X Purple Loading Dye to 5 $\mu$ L of reconditioned PCR pdt, run on either 2% agarose or 4-20% TBE acrylamide (to confirm successful elimination of heteroduplexes) at 180V for 40 minutes.

 DURATION

00:40:00