

# Segregation of plasmids in bacterial colonies

Tamara Matute, Isaac Núñez

## Abstract

**Citation:** Tamara Matute, Isaac Núñez Segregation of plasmids in bacterial colonies. **protocols.io**

dx.doi.org/10.17504/protocols.io.jbncime

**Published:** 10 Aug 2017

## Protocol

### Step 1.

Co-transform *E. coli* bacteria with two segregating plasmids (e.g. SEG10 and SEG11) and plate them on LB agar plates supplemented with arabinose and the antibiotics specific for each plasmid (e.g. carbenicillin + tetracycline).

#### 📌 NOTES

**Tamara Matute** 09 Aug 2017

Concentrations used:

Arabinose 10 mM

Chloramphenicol 10 µg/ml

Carbenicillin 100 µg/ml

Tetracycline 10 µg/ml

#### ■ ANNOTATIONS

**Fernan Federici** 13 Sep 2017

We use this protocol for competent cells:

([http://openwetware.org/wiki/TOP10\\_chemically\\_competent\\_cells](http://openwetware.org/wiki/TOP10_chemically_competent_cells))

### Step 2.

Incubate at 37°C overnight.

#### 📌 NOTES

**Tamara Matute** 09 Aug 2017

It is recommended to avoid long periods of growth because it tends to reduce the number of bacteria carrying the two plasmids within the colony (specially when carbenicillin is used).

### **Step 3.**

Check colonies by fluorescence imaging. The homogenous expression of both fluorescent proteins within the colony is a good indicator of co-transformed colonies.

### **Step 4.**

Pick a co-transformed colony and incubate these cells for 5 hours at 37°C with shaking (250 rpm) in liquid LB media supplemented with arabinose and the same specific antibiotics of each plasmid used in step 1 (e.g. carbenicillin + tetracycline).


### **Step 5.**

Harvest cells by centrifugation at 3000 rpm for 10 minutes and discard the supernatant.


### **Step 6.**

Re-suspend the cells in 5 ml of fresh LB media supplemented with the antibiotic in common to both SEG plasmids (e.g. chloramphenicol for SEG10 and SEG11).

If you wish to favor the formation of multipartite colonies go to step 8




If you wish to favor the formation of bipartite colonies go to optional step 7



(Optional Step) For bipartite colonies

**Step 7.**

Grow cells on liquid culture for extra 3h before going to step 8.



### Step 8.

Dilute this culture 1:5000 with fresh LB media supplemented with the antibiotic in common to both SEG plasmids (e.g. chloramphenicol for SEG10 and SEG11).

### Step 9.

Spread 10µl of this dilution on LB agar plates supplemented with the common antibiotic to both SEG plasmids (e.g. chloramphenicol for SEG10 and SEG11).

### NOTES

**Tamara Matute** 10 Aug 2017

Using 10µl produces ~50 colonies per plate. To vary this amount modify the dilution or the spread

volume at will. Fine tuning might be needed when SEG plasmids carry genetic components that impose high metabolic burden on the cell.

### **Step 10.**

Incubate at 37°C until colonies can be seen (8-12 hours).