

Removal of genomic DNA from RNA preparations (Thermo Scientific)

Maximilian Dietsch

Abstract

Removal of genomic DNA from RNA preparations

Citation: Maximilian Dietsch Removal of genomic DNA from RNA preparations (Thermo Scientific). **protocols.io**
dx.doi.org/10.17504/protocols.io.kzwcx7e

Published: 27 Nov 2017

Protocol

Step 1.

Add to an RNase free tube:

RNA	1 µg
10X reaction buffer with MgCl ₂	1 µl
DNase I, RNase-free	1 µl (1U)
Water	to 10 µl

Step 2.

Incubate at 37 °C for 30 min

🌡 **TEMPERATURE**

37 °C Additional info:

Step 3.

Add EDTA, Water and PCI and vortex thoroughly.

📦 **AMOUNT**

1 µl Additional info: EDTA

📦 **AMOUNT**

80 µl Additional info: Water

📦 **AMOUNT**

100 µl Additional info: PCI (phenol chloroform isoamyl alcohol)

Step 4.

Centrifuge for 10 min at 10.000 rpm and 4 °C

Step 5.

transfer the upper phase into a fresh tube and add 3 volumen EtOH/ 3M Natrumacetat (30:1, ph 5.2)

Step 6.

precipitate RNA over night at -20 °C