

# RNA extraction for the Betta splendens genome

Kailong Ma

## Abstract

This protocol is used to clarify the process of RNA extraction for our Betta splendens genome.

**Citation:** Kailong Ma RNA extraction for the Betta splendens genome. **protocols.io**

[dx.doi.org/10.17504/protocols.io.qvfdw3n](https://dx.doi.org/10.17504/protocols.io.qvfdw3n)

**Published:** 11 Jun 2018

## Protocol

### Sample preparation

#### Step 1.

1. Pour 1.5ml TRIZOL reagent.
2. For tissue samples, grind about 60mg with liquid nitrogen into powder and transfer the powder samples into the 2 ml tube contain of 1.5ml Trizol reagent.



#### REAGENTS

TRIZOL reagent 15596-026 by  
Invitrogen - Thermo Fisher

### Tissues lysis

#### Step 2.

Homogenize 2 minutes and place the sample at rest horizontally for 5 minutes to permit the complete dissociation of nucleoprotein complexes.

### Phase separation

#### Step 3.

1. Centrifuge at 12000×g for 5 minutes at 4°C. Transfer the supernatant to a new 2.0ml tube, add 0.3 ml of Chloroform / isoamyl alcohol(24:1) per 1.5 ml of Trizol reagent. Shake the tubes vigorously for 15 seconds.
2. Centrifuge at 12000×g for 10 minutes at 4°C. After centrifugation, the mixture should separates into three layers: the lower phenol-chloroform phase, an interphase, and an upper aqueous phase. RNA remains in the aqueous phase.



#### REAGENTS

Chloroform / isoamyl  
alcohol(24:1) 319988/W205702 by  
Sigma

## RNA precipitation

### Step 4.

1. Transfer the aqueous phase to a new 1.5mL tube; add equal volume of supernatant of isopropyl alcohol. Mixing well and place at -20°C for 2 hours for precipitation.
2. Centrifuge at 13600rpm for 20 minutes at 4°C and remove the supernatant.



#### REAGENTS

isopropyl alcohol w292907 by  
Sigma

## RNA washing

### Step 5.

1. Wash the RNA pellet with 1 ml 75% ethanol. Re-suspend the pellet and centrifuge at 13600rpm for 3 minutes at 4°C. Repeat this step again, Completely remove the ethanol without disturbing the pellet.
2. Air-dry the RNA pellet in the biosafety cabinet.

## Dissolve RNA

### Step 6.

Add 50µL of DEPC-treated water to dissolve the RNA pellet.



#### REAGENTS

DEPC-treated water AM9915G by  
Ambion