

Electroporation Protocol

New England Biolabs

Abstract

This protocol may be used with electrocompetent cells prepared by you according to [this protocol](#).

Citation: New England Biolabs Electroporation Protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.cruv6v

Published: 03 Feb 2015

Guidelines

Appropriate Antibiotics for Your Application

Antibiotics for Plasmid selection

Antibiotic	Working Concentration
Ampicillin	100 µg/ml
Carbenicillin	100 µg/ml
Chloramphenicol	33 µg/ml
Kanamycin	30 µg/ml
Streptomycin	25 µg/ml
Tetracycline	15 µg/ml

Electroporation Protocol

The electroporation protocol will vary depending on the strain so this protocol may need to be optimized. For control electroporation dilute pUC19 to 10 pg/µl with Milli-Q water.

Calculation:

If the culture was diluted 1000-fold when plated, the total cfu per ml is 1000 times the number of colonies counted. The cfu is divided by the amount of pUC19 (10 pg per ml)

$$\text{cfu/} \mu\text{g} = (\text{colonies counted} \times 1000) / (0.00001 \mu\text{g pUC19})$$

Before start

For control electroporation dilute pUC19 to 10 pg/µl with Milli-Q water.

Protocol

Step 1.

Turn on electroporator and set to 1.7-2.5 kv (optimize for strain), 200 ohms and 25 μ F

Step 2.

Place recovery SOC in 37°C water bath



. [SOC Media](#)

CONTACT: [New England Biolabs](#)

Step 2.1.

SOB Media



. [SOB Media](#)

CONTACT: [New England Biolabs](#)

Step 1.1.

2% tryptone

Step 1.2.

0.5% yeast extract

Step 1.3.

10 mM NaCl

Step 1.4.

2.5 mM KCl

Step 1.5.

10 mM MgCl₂

Step 1.6.

10 mM MgSO₄

Step 2.2.

20 mM glucose

Step 3.

Pre-warm LB-antibiotic plates at 37°C

Step 4.

Thaw cells on ice for 10 min or use freshly made cells



00:10:00

Step 5.

Place appropriate number of microcentrifuge tubes and 1 mm-electroporation cuvettes on ice

Step 6.

Flick the tube containing cells a few times to mix and add **25 μ l** to the microcentrifuge tubes



25 μ l Additional info:

Step 7.

Add **1 µl** of a 10 pg/µl DNA solution (in DI water) to the cells in the microcentrifuge tube

📄 **AMOUNT**

1 µl Additional info:

Step 8.

Transfer the DNA-cell mixture to the cold cuvette, tap on countertop 2X, wipe water from exterior of cuvette and place in the electroporation module and press pulse (**don't hold the button down**)

Step 9.

Immediately add **975 µl** of 37°C SOC, mix by pipetting up and down once and transfer to a 15 ml-falcon tube

📄 **AMOUNT**

975 µl Additional info:

Step 10.

Rotate in the 37°C incubator for 1 h

🕒 **DURATION**

01:00:00

Step 11.

Make appropriate dilutions

📌 **NOTES**

New England Biolabs 03 Feb 2015

When using 10 pg of DNA, make two dilutions:

Dilute 10 µl cells into 990 µl SOC and plate 100 µl. (1000-fold dilution)

Dilute 100 µl cells into 900 µl SOC and plate 100 µl. (100-fold dilution)

Step 12.

Incubate overnight at 37°C

🕒 **DURATION**

15:00:00

Warnings

The electroporation protocol will vary depending on the strain so this protocol may need to be optimized.