

Cooking Larvae Diet 👄

PLOS Genetics

Wen Aw1

¹z3314717@unsw.edu.au

dx.doi.org/10.17504/protocols.io.r6xd9fn



Wen Aw



EXT ERNAL LINK

https://doi.org/10.1371/journal.pgen.1007735

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Aw WC, Towarnicki SG, Melvin RG, Youngson NA, Garvin MR, Hu Y, Nielsen S, Thomas T, Pickford R, Bustamante S, Vila-Sanjurjo A, Smyth GK, Ballard JWO (2018) Genotype to phenotype: Diet-by-mitochondrial DNA haplotype interactions drive metabolic flexibility and organismal fitness. PLoS Genet 14(11): e1007735. doi: 10.1371/journal.pgen.1007735

PROTOCOL STATUS

Working

1	Weigh out agar 1 g per 100 ml of food to be prepared
2	Weigh out 0.022 g sugar, 2.719 g treacle, 7.938 g yeast and 7.321 g semolina per 100 ml of food for the 1:2 P:C diet, or 0.934 g sugar, 14.596 g treacle, 1.285 g yeast and 1.185 g semolina per 100 ml of food for the 1:16 P:C diet
3	Measure out water into large beaker, add small amount from large beaker to separate beaker
4	Add semolina to the small beaker and stir to ensure semolina is hydrated, leave aside
5	Add agar and sugar to large beaker and mix
6	Mark the volume on the beaker in pen
7	Place in microwave and bring to boil
8	Add treacle to large beaker, mix, and bring to boil

Add yeast to large beaker, mix, and bring to boil, continue heating until mixture stops over-boiling

10	Once the food will not boil over, add water to marked level to maintain the correct volume	
11	Boil for 10 minutes, mixing every 3 minutes, adding water to maintain volume	
12	Remove food from microwave and allow mixture to cool to 70°C	
13	Add semolina	
14	Allow food to cool to 60°C	
15	Add 1ml of 10% propionic acid, 1 ml of 1% phosphoric acid, and 0.4 ml of 25% nipagin dissolved in 100% ethanol per 100ml of food cooked	
16	Stir until semolina maintains suspension in the food	
17	Pour 25 ml of food per bottle	
18	Ensure semolina maintains suspension in the food mixture, stirring the food in the bottle if necessary	
19	Leave food to set and cool for 2 h	
20	Store food for no more than 1 week	
This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited		