Making your own electrocompetent cells

New England Biolabs

Abstract

This is the protocol for making two 250 ml cultures of electrocompetent cells

Citation: New England Biolabs Making your own electrocompetent cells. protocols.io

dx.doi.org/10.17504/protocols.io.crtv6m

Published: 05 Feb 2015

Before start

Sterile 10% glycerol (can be autoclaved) is needed for the washes. The volume of 10% glycerol needed is 2X the culture volume (for example, a 500 ml culture requires 1L of 10% glycerol).

Protocol

Step 1.

Innoculate 1 colony from a fresh plate of the strain to be made electrocompetent into **10 ml** of SOB in a 125 ml flask

■ AMOUNT

10 ml Additional info:

PROTOCOL

. SOB Media

CONTACT: New England Biolabs

Step 1.1.

2% tryptone

Step 1.2.

0.5% yeast extract

Step 1.3.

10 mM NaCl

Step 1.4.

2.5 mM KCl

Step 1.5.

10 mM MgCl2

Step 1.6.

10 mM MgSO4

Step 2.

Incubate for 16-18 hours at 37°C and 250 rpm

1

© DURATION

16:00:00

Step 3.

Pre-warm 2, 1 L flasks containing 250 ml each of SOB

■ AMOUNT

250 ml Additional info:

Step 4.

Add two drops of the overnight culture to each of the pre-warmed 250 ml flasks

■ AMOUNT

250 ml Additional info:

Step 5.

Shake at 37°C and 250 rpm until the cultures reach an OD600 of 0.5-0.7

■ AMOUNT

250 ml Additional info:

NOTES

New England Biolabs 03 Feb 2015

Be sure to turn on centrifuge and cool rotor to 4°C well in advance of harvesting cells. Be sure to place 1 L of 10% glycerol on ice well in advance of harvesting cells

Step 6.

Turn on centrifuge and cool rotor to 4°C well in advance of harvesting cells

Step 7.

Place **1** L of 10% glycerol on ice well in advance of harvesting cells

Step 8.

Place cultures on ice for 15 minutes. From this point on the cultures must be kept ice cold

O DURATION

00:15:00

Step 9.

Pour each 250 ml culture into chilled 500 ml (or 1000 ml) centrifuge bottles

Step 10.

Centrifuge at 5000 rpm for 10 min

O DURATION

00:10:00

Step 11.

Pour off the supernatant and aspirate any residual broth

Step 12.

Add **250 ml** of glycerol to each of the centrifuge bottles and completely suspend the cells by pipetting up and down

O DURATION

00:10:00

Step 13.

Centrifuge at 5000 rpm for 10 min

© DURATION

00:10:00

NOTES

New England Biolabs 03 Feb 2015

It is not necessary to aspirate

Step 14.

Pour off the supernatant, it is not necessary to aspirate

Step 15.

Add **250 ml** of glycerol to each of the centrifuge bottles and completely suspend the cells by pipetting up and down

Step 16.

Pour off the supernatant and suspend the cells in the **residual** glycerol by pipetting up and down

Step 17.

To freeze, add $100 \mu l$ of the culture to microcentrifuge tubes on ice

■ AMOUNT

100 µl Additional info:

NOTES

New England Biolabs 03 Feb 2015

At this point you can <u>electroporate</u> or freeze the cells away.

Step 18.

Once you have used all of the culture, transfer the tubes to dry ice for 10 minutes

O DURATION

00:10:00

Step 19.

Once the cultures are frozen, transfer them to a -80°C freezer