

Ligation of Sticky Ends with T4

Joshua Timmons

Abstract

Adapted T4 Protocol to reflect the one used by Northeastern_Boston for ligation of sticky end DNA.

Citation: Joshua Timmons Ligation of Sticky Ends with T4. **protocols.io**

dx.doi.org/10.17504/protocols.io.dqd5s5

Published: 05 Sep 2015

Protocol

Step 1.

Thaw the T4 DNA Ligase Buffer and resuspended at room temperature.

Step 2.

Set up the following reaction mixture at room temp:

10x T4 DNA Ligase Buffer	2 ul
Vector DNA	0.020 pmol
Insert DNA	0.060 pmol
Nuclease-free Water	to 20 ul
T4 DNA Ligase	1 ul

Step 3.

Gently mix the reaction by pipetting up and down and microfuge briefly.

Step 4.

For sticky ends, incubate at 22 C for 10 minutes.

Step 5.

Heat inactivate at 65 degrees C for 10 minutes.

 **DURATION**

00:10:00

Step 6.

Chill on ice and transform 1-5 µl of the reaction into 50 µl competent cells.