

Mar 21,  
2019

Working

## PBMC Isolation

Version 1

PLOS One

Dannielle Moore<sup>1</sup>

<sup>1</sup>[1] DST-NRF Centre of Excellence for Biomedical Tuberculosis Research; [2] South African Medical Research Council Centre for Tuberculosis Research; [3] Division of Molecular Biology and Human Genetics, Faculty of Medicine and Health Sciences, Stellenbosch University, Cape Town

[dx.doi.org/10.17504/protocols.io.yfvftn6](https://doi.org/10.17504/protocols.io.yfvftn6)

Dannielle Moore

[1] DST-NRF Centre of Excellence for Biomedical Tuberculosis...



## ABSTRACT

Isolation of human peripheral blood mononuclear cells (PBMCs) using the Ficoll density gradient method.

## EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0213832>

## PROTOCOL STATUS

**Working**

## GUIDELINES

Unless otherwise stated, sample processing was performed in a biosafety cabinet under sterile conditions.

## MATERIALS

NAME	CATALOG #	VENDOR
Falcon Tube (50 mL)		Fischer Scientific
Disposable pasteur pipettes	EA61.1	Carl Roth
Ficoll-Paque PLUS density gradient media	17144002	Ge Healthcare
1x Phosphate-Buffered Saline	04-479Q	Lonza
Fetal Bovine Serum	SH30088.02	HyClone
Trypan Blue Solution 0.4% Sterile-filtered	T8154	Sigma Aldrich
Vac 9mL Sodium Heparin Green	VGRV455051	Lasec

## SAFETY WARNINGS

## BEFORE STARTING

Collect human whole peripheral blood in Sodium Heparin tubes – This step is done by a professional healthcare worker in a medical examination room.

- 1 In a 50mL Falcon tube, add 15mL of Ficoll-Histopaque Plus media
- 2 In a separate 50mL Flacon tube, dilute peripheral blood in 1:1 ratio with 1x phosphate-buffered saline to a max volume of 35mL (max starting volume of blood is 18mL)

Note: if a larger volume of blood is required - the blood must be split across separate falcon tubes and isolated PBMCs combined following isolation procedure, prior to cell counting

- 3 Gently layer the diluted blood from step 3 onto the Ficoll from step 2. This can be done free-hand or using a graduated pipette and pipette man
- 4 Once the blood has been layered and tubes secured, remove from hood and insert into benchtop centrifuge.
- 5 Spin at 400xg for 25 min at room temperature with the accelerator and brake OFF.
- 6 Once centrifuge has stopped, remove tubes and place in hood for further processing.
- 7 Use a sterile Pasteur pipette, carefully remove the upper plasma layer until 5cm above opaque PBMC band. In a circular motion, collect the PBMC band at the Ficoll interface and transfer into a new 50mL Falcon tube

Note: if a large volume of blood was processed and blood split into several tubes, collect and decant all PBMC bands into single falcon tube.

- 8 Wash isolated PBMCs twice in 50mL of PBS. Centrifuge cell suspension at 400xg for 10 min at room temperature with the brake and accelerator set to max.
- 9 Dilute cells in 1:10 ratio with PBS and trypan blue. Count the cells using haemocytometer and microscope. Record observed cell number and cell viability



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited