



## Single-cell mapping of lineage and identity via CellTagging

Samantha Morris<sup>1</sup>, Brent A. Biddy<sup>1</sup>, Wenjun Kong<sup>1</sup>, Kenji Kamimoto<sup>1</sup>, Chuner Guo<sup>1</sup>, Sarah Wayne<sup>1</sup>, Tao Sun<sup>1</sup>

<sup>1</sup>Washington University, Saint Louis

[dx.doi.org/10.17504/protocols.io.u9yez7w](https://doi.org/10.17504/protocols.io.u9yez7w)

Samantha Morris

### ABSTRACT

Single-cell technologies are offering unprecedented insight into complex biology, revealing the behavior of rare cell populations that are typically masked in bulk population analyses. One current limitation of single-cell approaches is that lineage relationships are lost as a result of cell processing, restricting interpretations of the data collected. Elegant computational approaches have been developed in an effort to infer these missing observations, but it remains a challenge to reconstruct true reprogramming trajectories using these tools. Although sophisticated lineage tracing solutions to connect cell history with fate are emerging, these protocols are either not compatible with high-throughput scRNA-seq, or require genome editing strategies that are not readily deployed in some systems. Here, our protocol describes a single-cell resolution clonal tracking approach, 'CellTagging', based on combinatorial cell indexing, permitting the parallel capture of lineage information and cell identity. CellTagging integrates with high-throughput single-cell RNA-sequencing, where iterative rounds of cell labeling enable the construction of multi-level lineage trees. This straightforward lentiviral-labeling approach can be applied to an array of cell biological applications to simultaneously profile lineage and identity, at single-cell resolution.

### TAGS

Single Nuclei RNA-Seq

Reprogramming

Show tags

### THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

<https://www.biorxiv.org/content/early/2017/04/28/127860>

### PROTOCOL STATUS

#### In development

We are still developing and optimizing this protocol

### MATERIALS

NAME	CATALOG #	VENDOR
Stellar Competent Cells	636763	Takarabio
CellTag pooled library V1	115643	addgene
CellTag pooled library V2	115644	addgene
CellTag pooled library V3	115645	addgene

### SAFETY WARNINGS

For generation of lentivirus, the appropriate institutional regulations must be satisfied.

#### Amplification of CellTag library

1

Thaw Stellar Competent Cells in an ice bath just before use.

- 2 After thawing, mix 100 uL of cells with 10-50 ng of CellTag V1 DNA in a 1.5-mL microcentrifuge tube.
- 3 Place transformation mixture on ice for 30 minutes.
- 4 Heat shock the cells for exactly 60 seconds at **42 °C**.
- 5 Place tube on ice for 1-2 minutes.
- 6 Add SOC medium to bring the final volume to 1000 uL.
- 7 Incubate by shaking (~250 rpm) for 1 hour at **37 °C**.
- 8 Take 6 uL of the recovery. Use 1 uL of this to prepare a 1:10 and then a 1:100 dilution. Plate both dilutions and the rest of the 5 uL of recovery onto three separate LB/Amp plates. Grow overnight at **37 °C**.
- 9 Add the rest of the recovery to 500 mL of LB+Amp. Grow overnight while shaking (~250 rpm) at **37 °C**.
- 10 After counting the CFUs from the LB/Amp plates the next day, calculate the transformation efficiency.
- 11 Harvest the cells from the liquid culture and use multiple Maxiprep columns (Sigma or Zymogen) to purify the library.
- 12



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited