

Digestion BclI

Antoine ZILLER

Abstract

Citation: Antoine ZILLER Digestion BclI. **protocols.io**

dx.doi.org/10.17504/protocols.io.ecvbaw6

Published: 05 Jan 2016

Protocol

Step 1.

Add:

- nuclease-free water qsp 50 μL > 43 μL
- 10X Buffer G 5 μL
- 1 μg DNA (1 $\mu\text{g}/\mu\text{L}$)
- 1 μL BclI 10U/ μL > 10U/1 μg DNA

Mix gently and spin down for a few seconds.

****Star Activity :** An excess of BclI (20 U/ μg DNA x 1 hour) may result in star activity

Step 2.

Incubate at 55°C for 1-16 hours**.

It's possible to use extended digestion for 16 hours by adding 0.5U for 1 μg of DNA in 50 μL

 **DURATION**

01:30:00

Step 3.

- Inactivated at 80°C in 20 min.

> For only 10 U of enzyme

or - adding 0.5 M EDTA, pH 8.0 (#R1021), to achieve a 20 mM final concentration

Step 4.

1. Measure the volume of the DNA sample.

> 50 μL

2. Add 1/10 volume of sodium acetate, pH 5.2, (final concentration of 0.3 M) - These amounts assume that the DNA is in TE only; if DNA is in a solution containing salt, adjust salt accordingly to achieve the correct final concentration.

> 5 μL

3. Mix well.

4. Add 2 to 2.5 volumes of cold 100% ethanol (calculated after salt addition).
> 110 (2V)
5. Mix well.
6. Place on ice or at -20 degrees C for >20 minutes.
7. Spin a maximum speed in a microfuge 10-15 min.
8. Carefully decant supernatant.
9. Add 1 ml 70% ethanol. Mix. Spin briefly. Carefully decant supernatant.
10. Air dry or briefly vacuum dry pellet.
11. Resuspend pellet in the appropriate volume of TE or water.