

Low Gel Temperature Agarose (LGTA) Media

G Jason Smith, April Woods

Abstract

Use for solid phase culture and chemical selection of *Pseudo-nitzschia spp* and other recalcitrant diatoms

Citation: G Jason Smith, April Woods Low Gel Temperature Agarose (LGTA) Media. **protocols.io**

dx.doi.org/10.17504/protocols.io.gj4buqw

Published: 30 Nov 2016

Protocol

REAGENTS:

Step 1.

- 0.2 µm filtered seawater (FSW)
- L1 Medium Kit (ncma.bigelow.org)
- SeaPrep™ Agarose (Lonza Cat#: 50302, 25 gm)

PREPARATION: 1% (w/v) LGTA

Step 2.

- a. Add 100 mL FSW to autoclave compatible screw top bottle with stir bar
- b. To stirring FSW add:
 - 1.00 gm SeaPrep agarose to stirring FSW, *avoid clumping*
 - 0.10 mL L1 - NaH₂PO₄ stock (36.2 µM final, 1X L1 conc)
 - 0.05 mL L1- NaNO₃ stock (441 µM final, 0.5X L1 conc)
 - 0.20 mL L1 - Na₂SiO₃ stock (212 µM final, 2X L1 conc)
 - 0.10 mL L1 - Trace Element stock (1X L1 conc)
- c. Autoclave mixture to fully melt agarose.
- d. Cool in water bath to 20° to 25°C with occasional stirring
- e. Add 0.05 mL L1 - Vitamin stock (standard conc.), gently mix

NOTE: LGTA stocks can be held in molten form at ≥ 20°C or allowed to solidify at ≤ 15°C

and remelted in microwave prior to cooling for use. *KEEP STERILE*

PLATING CULTURES:

Step 3.

- a. Aliquot cooled, molten LGTA into sterile culture tubes (e.g. Falcon 2059).
 - Volume should be sufficient to form thin layer (<5mm) in target plates or wells to enable in-gel imaging.
 - *Recommend 0.2 mL per 24 well plate well, 2.5 mL for 30mm diameter petri plates, and 5 mL for 60 mm diameter plates.*
- b. Add selection reagents as needed.
- c. Add cell suspension **into** LGTA.
- d. Gently vortex to disperse cells.
- e. Transfer to culture plate avoiding introducing bubbles.
- f. Cool and incubate at 15°C to set agarose.
- g. Transfer to standard growth conditions as needed
 - For species with high growth temperatures solidified LGTA can be transferred to appropriate incubation temperature after solidifying.
 - For long term solid phase culture, the LGTA gel can be overlaid with L1 media.

NOTE: Good gel strength is obtained at $\geq 0.75\%$ (w/v) LGTA and optical clarity obtained at $\leq 2\%$ LGTA).

We have had good luck with *Pseudo-nitzschia* growing in $\geq 0.75\%$ to $\leq 1\%$ LGTA.

Concentration of stock LGTA can be adjusted to accommodate larger culture volumes but a

1:1 (v/v) mixture of culture to LGTA should be considered the maximum mixing ratio for good gel formation.

