Preparation of Cas9 Protein

Ashwin Banfal, Brett T. Staahl, Steve Lin

Abstract

This is a protocol from the <u>Doudna Lab</u> for preparation of Cas9 protein.

Citation: Ashwin Banfal, Brett T. Staahl, Steve Lin Preparation of Cas9 Protein. protocols.io

dx.doi.org/10.17504/protocols.io.rczd2x6

Published: 21 Aug 2018

Guidelines

The protocol workflow is as follows:

1. Expression of Cas9 in Rosetta2 E. coli

- Transformation (<u>Steps 1-5</u>)
- Starter culture (Steps 6-8)
- Expression (<u>Steps 9-21</u>)

2. Purification of Cas9

- Sonication (Steps 22-29)
- HisTrap (Steps 30-40)
- TEV (<u>Steps 41-45</u>)
- Heparin (<u>Steps 46-53</u>)
- BPTrap/Superdex 200 (<u>Steps 54-65</u>)

Before start

Materials

TY media 1 L
Tryptone 10 g
Yeast extract 6 g
NaCl 10 g

Dissolve in 400 mL

Rinse flasks thoroughly with MQ water Add 600 mL MQ water, media

Autoclave (45', no dry time)

Buffers

Lysis	1L, pH7.5
20mM HEPES	20 ml of 1M
1M KCI	74.55 g
10% glycerol	100 ml
5mM TCEP	1.43 g
10mM imidazole	680 mg

Elution	500mL, pH7.5
20mM HEPES	10 ml of 1M
0.1M KCI	50 ml of 1M
10% glycerol	50 ml
5mM TCEP	717 mg
300mM imidazole	10.2 g

Ion Exchange
buffer A1L, pH7.5 (low
salt)20mM HEPES20 ml of 1M300mM KCl300 ml of 1M

10% glycerol 100 ml 1mM TCEP 287 mg Ion Exchange 1L, pH7.5 (high buffer B salt)

 20mM HEPES
 20 ml of 1M

 1M KCl
 74.55 g

 10% glycerol
 100 ml

 1mM TCEP
 287 ml

Gel Filtration1L, pH7.520mM HEPES20 ml of 1M150mM KCl150 ml of 1M10% glycerol100 ml1mM TCEP287 mg

Materials

Please see before starting in Guidelines section for materials. by Contributed by users

Protocol

Expression of Cas9 in Rosetta2 E. coli - Transformation

Step 1.

Mix 1 μl of 10 ng/uL DNA + 10 μl Rosetta 2 cells (pRAR-rare codons).

■ AMOUNT

1 μl: 10 ng/uL DNA

■ AMOUNT

10 μl : Rosetta 2 cells

Expression of Cas9 in Rosetta2 E. coli - Transformation

Step 2.

Heat shock: 10' ice, 45" 42°C, 5' ice.

© DURATION 00:10:00 : ice © DURATION 00:00:45 : 42°C © DURATION 00:05:00 : ice

Expression of Cas9 in Rosetta2 E. coli - Transformation

Step 3.

Recover: Add 300 µl LB, 50' 37°C on shaker.

■ AMOUNT 300 µl : LB

↓ TEMPERATURE
37 °C : Shaker
○ DURATION
00:50:00 : Shaker

Expression of Cas9 in Rosetta2 E. coli - Transformation

Step 4.

Spin down (6000 xg, 1'), and remove 150 µl LB.

O DURATION

00:01:00 : Spin down

Expression of Cas9 in Rosetta2 E. coli - Transformation

Step 5.

Plate 100 µl, 30 µl on 2 AMP plates.

Expression of Cas9 in Rosetta2 E. coli - Starter culture

Step 6.

Pick 6 colonies to add to 20 ml LB + Amp (Cf=100 ug/ml) starter culture.

AMOUNT

20 ml : LB + Amp (Cf=100 ug/mL) starter culture

Expression of Cas9 in Rosetta2 E. coli - Starter culture

Step 7.

Grow overnight.

Expression of Cas9 in Rosetta2 E. coli - Starter culture

Step 8.

Next day: Store XL1 blue 20% glycerol stocks in -80°C.

↓ TEMPERATURE

-80 °C : Storage

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 9.

Add AMP to the flasks (Cf=100 ug/ml).

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 10.

Add 20 ml starter culture in the morning,

■ AMOUNT

20 ml: starter culture

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 11.

Grow to OD600=0.6 (150 rpm, 4 hrs, 37°C).

TEMPERATURE

37 °C : Growing

O DURATION

04:00:00 : Growing

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 12.

Set shaker temp to 16°C.

■ TEMPERATURE

16 °C : Shaker temperature

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 13.

*Save 1 ml for SDS-PAGE: uninduced.

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 14.

Add 0.5 mM IPTG.

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 15.

Grow overnight (150 rpm, at least 16 hrs, 16°C).

↓ TEMPERATURE
 16 °C : Growing
 ○ DURATION
 16:00:00 : Growing

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 16.

Harvest cells before lunch: Spin down (4000 rpm, 20', 4°C).

↓ TEMPERATURE
4 °C : Spin down
⑤ DURATION

00:20:00 : Spin down

NOTES

*Keep the cells cold from here.

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 17.

Completely resuspend pellets in 20 ml Lysis Buffer + 1mM PMSF per bottle.

■ AMOUNT

20 ml: Lysis Buffer + 1mM PMSF

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 18.

Use additional 20 ml Lysis Buffer to wash all bottles.

■ AMOUNT

20 ml: Lysis Buffer

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 19.

Transfer to 50 ml-tubes.

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 20.

*Save 20 µl for SDS-PAGE: cell pellet

Expression of Cas9 in Rosetta2 E. coli - Expression

Step 21.

Freeze the cells at -80°C.

↓ TEMPERATURE

-80 °C: Freezing cells

Purification of Cas9 - Sonication

Step 22.

Thaw cells in water.

Purification of Cas9 - Sonication

Step 23.

Push beaker with a stir bar to the bottom of ice bucket.

Purification of Cas9 - Sonication

Step 24.

Add 50 ml Lysis Buffer and 1 protease inhibitor tablet/50 ml cell. Dissolve the tablets.

■ AMOUNT

50 ml: Lysis Buffer

Purification of Cas9 - Sonication

Step 25.

Transfer the cells to the beaker.

Purification of Cas9 - Sonication

Step 26.

Sonicate while stirring: 8' total, 10" on, 10" off, level 4. Lysate should turn fluid and thin. Repeat 8 min cycle, if necessary.

© DURATION 00:08:00 : total © DURATION 00:00:10 : on/off

Purification of Cas9 - Sonication

Step 27.

Transfer to Oak Ridge tubes (30 ml/tube). *Balance within +/- 0.01g.

Purification of Cas9 - Sonication

Step 28.

Pellet cell debris (18,000 rpm, 30', 4°C). Prepare HisTrap purification.

HisTrap

CV=5mL, FR=2mL/min

Buffer A: Lysis buffer, Buffer B: Elution buffer

↓ TEMPERATURE 4 °C : Pellet cell debris

O DURATION

00:30:00 : Pellet cell debris

Purification of Cas9 - Sonication

Step 29.

*Save 20 µl supernatant for SDS-PAGE: cell lysate

Purification of Cas9 - HisTrap

Step 30.

Filter sample.



HisTrap

CV=5ml, FR=2ml/min

Buffer A: Lysis buffer, Buffer B: Elution buffer

Purification of Cas9 - HisTrap

Step 31.

Set-up Akta:

- 1. Turn on UV lamp.
- 2. Wash Akta with water: pump wash, and flow 50% A, 50% B.
- 3. Pump wash.
- 4. Flow 50% A, 50% B.
- 5. Flow 100% A, until steady baseline. Set baseline to zero.

Purification of Cas9 - HisTrap

Step 32.

Apply sample with pump, in the cold room:

- 1. Fill tubing with water, connect column.
- 2. Wash: 5CV water.
- 3. Equilibrate: 5CV lysis buffer.
- 4. Apply sample to the column. Collect flow-through: Ni-sample loading.

Re-apply flow-through to the column, if necessary.

Purification of Cas9 - HisTrap

Step 33.

Purification:

In the cold room:

First wash: 15CV lysis buffer. Collect flow-through: Ni-wash1.

Collect flow-through: Ni-wash1. *Save 20µl for SDS-PAGE.

Purification of Cas9 - HisTrap

Step 34.

Wash2: 3CV lysis buffer + 0.1% Triton X.

Collect flow-through: Ni-wash2.

*Save 20 µl for SDS-PAGE.

Purification of Cas9 - HisTrap

Step 35.

Purification using Akta:

First, attach column to Akta (FR=0.9 ml/min).

Purification of Cas9 - HisTrap

Step 36.

Wash: 5CV lysis buffer, until steady UV baseline.

Purification of Cas9 - HisTrap

Step 37.

Elute: 25CV elution buffer, 100% gradient.

Purification of Cas9 - HisTrap

Step 38.

Collect fractions (Fraction size=2.5 ml).

Purification of Cas9 - HisTrap

Step 39.

Clean-up:

- 1. Wash: 5CV water
- 2. Store: 5CV 20% ethanol
- 3. Detach column.
- 4. Wash Akta with water.

Purification of Cas9 - HisTrap

Step 40.

Pool appropriate fractions. Rinse collection tubes with elution buffer.

*Save 10 µl for SDS-PAGE: Ni-elution

Purification of Cas9 - TEV

Step 41.

Concentrate with 30K Amicon to 3 ml volume. Collect flow-though: Ni-Amicon FT.

Purification of Cas9 - TEV

Step 42.

Add 1.25 mg TEV/3 ml sample, and mix with pipet.

AMOUNT

1.25 mg: TEV

Purification of Cas9 - TEV

Step 43.

Incubate on shaker, at 4°C overnight.

↓ TEMPERATURE

4 °C: Incubation on shaker

© DURATION

16:00:00 : overnight incubation

Purification of Cas9 - TEV

Step 44.

Next day: *Save 10 µl for SDS-PAGE: **postTEV**

Purification of Cas9 - TEV

Step 45.

Add 25 ml of buffer A to sample.

■ AMOUNT

25 ml: Buffer A

Purification of Cas9 - Heparin

Step 46.

Heparin

CV=5 mL, FR=1 mL/min

Buffer A: IEX low salt buffer, Buffer B: IEX high salt buffer

Purification of Cas9 - Heparin

Step 47.

Set-up Akta (FR=2 ml/min).

Purification of Cas9 - Heparin

Step 48.

Apply sample with pump, in the cold room:

- 1. Equilibrate: 10CV buffer A.
- 2. Apply sample to the column. Collect flow-through: Hep-sample loading.

Purification of Cas9 - Heparin

Step 49.

Purification:

- 1. Attach column to Akta (FR=0.9 ml/min).
- 2. Wash: 15CV buffer A, until steady baseline.
- 3. Collect flow-through: **Hep-wash**. *Save 20 µl for SDS-PAGE.
- 4. Elute: 20CV buffer B, 100% gradient.

5. Collect fractions (Fraction size=2 ml).

Purification of Cas9 - Heparin

Step 50.

Clean-up (FR=2ml/min):

- 1. Wash: buffer B, until steady conductivity baseline for over 30 ml.
- 2. Store: 5CV 20% ethanol.
- 3. Detach column.
- 4. Wash Akta with water.

Purification of Cas9 - Heparin

Step 51.

Pool appropriate fractions. Rinse collection tubes with buffer B.

*Save 10 µl for SDS-PAGE: **Hep-elution**

Purification of Cas9 - Heparin

Step 52.

Concentrate sample to 1 ml. Collect flow-through: Hep-Amicon FT.

Purification of Cas9 - Heparin

Step 53.

Leave concentrator on ice, at 4°C overnight.

■ TEMPERATURE

4 °C: Concentrator on ice

Purification of Cas9 - MBPTrap/Superdex 200

Step 54.

MBPTrap/Superdex 200

CV=120 mL, FR=0.5-1 mL/min

Previous day:

- 1. Set-up Akta (FR=2 mL/min):
- 2. Prepare columns:
- a. MBPTrap wash: 5CV water, 5CV gel filtration buffer (FR=2 mL/min), in the cold room.
- b. Attach MBPTrap and Superdex200 to Akta.
- c. Equilibrate overnight: 2CV gel filtration buffer (FR=0.3 mL/min).
- *Set end timer volume=240 mL, max pressure alarm=0.25 MPa.

Purification of Cas9 - MBPTrap/Superdex 200

Step 55.

Concentrate sample to 500 µl. Rinse concentrator with 500 µl gel filtration buffer.

■ AMOUNT

500 μl : Gel filtration buffer

Purification of Cas9 - MBPTrap/Superdex 200

Step 56.

Spin down sample (11000 rpm, 4', 4°C). Transfer supernatant to a new 1.5 ml-tube.

TEMPERATURE
4 °C : Spin down
C DURATION

00:04:00 : Spin down

Purification of Cas9 - MBPTrap/Superdex 200

Step 57.

To inject sample, first, clean injection port, sample loop with water, gel filtration buffer.

Purification of Cas9 - MBPTrap/Superdex 200

Step 58.

Now, draw up sample with 1 ml-syringe, and replace the syringe on the injection port.

Purification of Cas9 - MBPTrap/Superdex 200

Step 59.

[LOAD] Inject sample to the loop with syringe.

Purification of Cas9 - MBPTrap/Superdex 200

Step 60.

Purification

- 1. [INJECT] Apply sample to the column (3 ml).
- 2. [LOAD] Flow 1CV gel filtration buffer. Cas9 elutes at 65-85 ml.
- 3. Collect fractions (Fraction size=2 ml).

Purification of Cas9 - MBPTrap/Superdex 200

Step 61.

Clean-up (FR=2 ml/min):

- 1. Detach MBPTrap.
- 2. Wash MBPTrap: 10CV 10mM maltose + gel filtration buffer, in the cold room.
- 3. Store MBPTrap: 5CV 20% ethanol.
- 4. Wash Superdex200: 2CV water (Long-term storage in 20% ethanol).
- 5. Detach Superdex200.
- 6. Wash Akta with water.

Purification of Cas9 - MBPTrap/Superdex 200

Step 62.

Pool appropriate fractions. Rinse collection tubes. *Save 5 µl for SDS-PAGE: Cas9

Purification of Cas9 - MBPTrap/Superdex 200

Step 63.

Concentrate to desired concentration (40 uM).

*Cas9

MW=158,441.4 g/mol

Molar absorption coefficient=120,575 L/mol·cm

Purification of Cas9 - MBPTrap/Superdex 200

Step 64.

Aliquot 20 µl each into 1.5 ml-tubes.

Purification of Cas9 - MBPTrap/Superdex 200

Step 65.

Snap-freeze in LN2, and store at -80°C.

↓ TEMPERATURE

-80 °C: Storage

Warnings

Please refer to the SDS (Safety Data Sheet) for safety warnings and hazard information.