

# Preparation of Virus DNA from Seawater for Metagenomics

Matthew Sullivan Lab

## Abstract

**Citation:** Matthew Sullivan Lab Preparation of Virus DNA from Seawater for Metagenomics. **protocols.io**

dx.doi.org/10.17504/protocols.io.c28yhv

**Published:** 04 Jan 2016

## Guidelines

Needed:

- 0.2 µm filtered seawater with iron chloride
- Mg-EDTA-Ascorbate (Oxalate) buffer
- DNase I
- EDTA
- EGTA
- Centrifugal concentrators
- Resuspension buffer
- Wizard prep column
- Wizard prep resin
- Quant-iT dsDNA Pico Green assay kit

## Materials

Quant-iT dsDNA Pico Green assay kit (Invitrogen) [P7589](#) by [Life Technologies](#)

Wizard Minicolumns [A7211](#) by [Promega](#)

Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-100 membrane [910024](#) by [Emd Millipore](#)

## Protocol

### Step 1.

Precipitate 0.2µm filtered seawater with iron chloride and store the filters at 4°C

#### ⊕ NOTES

**VERVE Team** 14 Jul 2015

See [Ferric Chloride Precipitation of Viruses from Seawater protocol](#) for more details

### Step 2.

Resuspend the precipitated virus in Mg-EDTA-Ascorbate (or Oxalate) buffer just prior to DNA extraction (*also see Ferric Chloride Resuspension Buffer protocol*)

### Step 3.

Treat the resuspended virus preparation with DNase I for 2 hours at room temperature

#### DURATION

02:00:00

### Step 4.

Inactivate the enzyme with EDTA and EGTA

#### NOTES

**VERVE Team** 16 Jun 2015

Previous work in the Tucson Marine Phage Lab has shown that it is not necessary to perform CsCl gradient centrifugation on the virus preparation prior to DNA extraction. The iron chloride precipitation followed by DNase I treatment removes the majority of bacterial DNA that could contaminate the final metagenomic sequence. CsCl may be used, but for deep or other samples where viruses are limiting, the recovery of DNA may be significantly diminished by performing this step.

**VERVE Team** 14 Jul 2015

See [DNase I Treatment Protocol](#) for more details.

### Step 5.

Concentrate the resuspended virus preparation using centrifugal concentrators down to 3-5 ml total volume if performing CsCl cleanup or to 1-2 ml total volume if proceeding directly to DNA extraction

#### REAGENTS

Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-100 membrane [910024](#) by [Emd Millipore](#)

#### NOTES

**VERVE Team** 16 Jun 2015

The Sullivan Lab recommends Amicon Ultra, 100kDa MWCO, either 15 ml or 5 ml capacity (Millipore UFC910024 or UFC810024), following manufacturer's directions.

### Step 6.

Remove the concentrate to a fresh tube

### Step 7.

Rinse the membrane with an additional 0.5-1 ml resuspension buffer by pipetting up and down along the membrane

#### PROTOCOL

#### . [Resuspension Buffer P1](#)

CONTACT: [Bonnie Poulos](#)

#### Step 7.1.

Dissolve 6.06g Tris Base in 800 mL MilliQ water

#### REAGENTS

Tris Base [BP152-1](#) by [Fisher Scientific](#)

#### ANNOTATIONS

**Chris Upton** 12 Oct 2015

Isn't this usually made from stock Tris and EDTA solutions?

**Bonnie Poulos** 12 Oct 2015

It can be made either way depending on what you have available in the lab. Some find it easier to use stock solutions, but it is not necessary as long as final pH is adjusted.

#### Protocol

##### Step 7.2.

Add 3.72g EDTA disodium salt, dihydrate to the 800 mL Tris base and stir to dissolve

 **AMOUNT**

6 g Additional info:

 **REAGENTS**

EDTA, disodium salt, dihydrate [S312-500](#) by [Fisher Scientific](#)

#### Protocol

##### Step 7.3.

Adjust the pH to 8.0 with HCl

#### Protocol

##### Step 7.4.

Adjust the volume to 1 liter with MilliQ water

#### Protocol

##### Step 7.5.

Add 100mg RNase A per liter of buffer P1

 **AMOUNT**

100 mg Additional info:

##### Step 8.

Pool with the concentrated sample.

##### Step 9.

At this point, either layer this onto premade CsCl gradients and isolate the virus fraction according to the protocol [Cesium Chloride Purification of Viruses](#) (or alternatively, [Cesium Chloride DNA Extraction of Viruses using Wizard Prep Columns](#)) or proceed directly to Wizard Prep purification of DNA.

##### Step 10.

Extract the DNA using Wizard Prep columns and resin

#### NOTES

**VERVE Team** 14 Jul 2015

See [Cesium Chloride DNA Extraction of Viruses using Wizard Columns protocol](#) for more details

##### Step 11.

Calculate the amount of DNA recovered using Quant-iT dsDNA Pico Green assay kit following manufacturer's directions

 **REAGENTS**

Quant-iT dsDNA Pico Green assay kit (Invitrogen) [P7589](#) by [Life Technologies](#)

#### NOTES

**VERVE Team** 26 Jun 2015

The community virus DNA is now ready for metagenomic library preparation or for linker amplification if more DNA is needed for sequencing (see Linker Amplification protocol and Bar-Coded Oligos for LA if needed).

**VERVE Team** 14 Jul 2015

See [Quant-iT Pico Green dsDNA Assay Protocol](#) for more details