# **Blunting Protocol for NEB PCR Cloning Kit (E1202)**

# **New England Biolabs**

# **Abstract**

This is the blunting protocol for NEB PCR Cloning Kit (E1202)

Citation: New England Biolabs Blunting Protocol for NEB PCR Cloning Kit (E1202). protocols.io

dx.doi.org/10.17504/protocols.io.crpv5m

Published: 03 Feb 2015

## **Guidelines**

Reaction volume may be scaled up or down as necessary.

PCR generated DNA must be purified before blunting by using a commercial purification kit, phenol extraction/ethanol precipitation, or gel electrophoresis.

Restriction enzyme digested DNA can be blunted directly without purification. The Blunt Enzyme Mix has been optimized in <u>Blunting Buffer</u>, but is also active in NEBuffers 1.1, 2.1, 3.1, and CutSmart™ Buffer in addition to NEBuffers 1-4, BamHI, EcoRI and DpnII unique buffers when supplemented with dNTPs and dithiothreitol. There is a small reduction in ligation fidelity in these buffers. Transformation efficiency is lowest in NEBuffer 1 and 1.1 where the total yield is about 50% of optimum.

## **Before start**

Reaction volume may be scaled up or down as necessary.

# **Materials**

NEB PCR Cloning Kit - 20 rxns <u>E1202S</u> by <u>New England Biolabs</u>

# **Protocol**

## Step 1.

Mix the following components in a sterile microfuge tube:

**PROTOCOL** 

. E1202 Blunting Mixture

**CONTACT: New England Biolabs** 

Step 1.1.

Purified DNA (up to 5 μg): 1-19 μl

Step 1.2.

10X Blunting Buffer: 2.5 μl

**■** AMOUNT

3 μl Additional info:

Step 1.3.

1 mM dNTP Mix: **2.5 μl** 

**■** AMOUNT

3 µl Additional info:

Step 1.4.

Blunt Enzyme Mix: 1.0 µl

**■** AMOUNT

1 µl Additional info:

Step 1.5.

Sterile dH20 up to 25 µl

# Step 2.

Reactions containing restriction enzyme digested DNA are incubated at room temperature for 15 minutes. Reactions with sheared/nebulized DNA or PCR products are incubated at room temperature for 30 minutes.

#### NOTES

## New England Biolabs 03 Feb 2015

PCR generated DNA must be purified before blunting by using a commercial purification kit, phenol extraction/ethanol precipitation, or gel electrophoresis.

## New England Biolabs 03 Feb 2015

Restriction enzyme digested DNA can be blunted directly without purification. The Blunt Enzyme Mix has been optimized in <u>Blunting Buffer</u>, but is also active in NEBuffers 1.1, 2.1, 3.1, and CutSmart™ Buffer in addition to NEBuffers 1-4, BamHI, EcoRI and DpnII unique buffers when supplemented with dNTPs and dithiothreitol. There is a small reduction in ligation fidelity in these buffers. Transformation efficiency is lowest in NEBuffer 1 and 1.1 where the total yield is about 50% of optimum.

## Step 3.

Immediately inactivate enzyme in the blunting reaction by heating at 70°C for 10 minutes.

**O DURATION** 

00:10:00

#### Step 4.

Proceed directly to the ligation step using the Quick Ligation Kit (NEB #M2200) or standard T4 DNA Ligase (NEB #M0202).

## NOTES

# New England Biolabs 03 Feb 2015

Blunt ligation reactions using standard T4 DNA Ligase should be incubated overnight at room temperature.