

# Genotyping the ADAMTS5-P (Pfizer) mouse strain

Bianca Hemmeryckx, Dries Bauters, H. Roger Lijnen

## Abstract

ADAMTS5, the main aggrecanase, is a protein that seems to play a functional role in the development of brown and white adipose tissue (WAT) and browning of WAT. These published observations were made using the ADAMTS5-P mice originally generated by Pfizer, in collaboration with Lexicon Genetics (The Woodlands, TX, USA).

**Citation:** Bianca Hemmeryckx, Dries Bauters, H. Roger Lijnen Genotyping the ADAMTS5-P (Pfizer) mouse strain. **protocols.io**

dx.doi.org/10.17504/protocols.io.kaucsew

**Published:** 14 Oct 2017

## Guidelines

-Tails digested in lysis buffer containing proteinase K can be stored at -20°C before continuing the protocol.

## Protocol

### Step 1.

Cut a very small piece of the ear

### Step 2.

Collect them in an labeled eppendorf tube filled with 800 µl lysis buffer supplemented with 8µl 10 mg/ml proteinase K solution

\* Lysis buffer (1L):

-40 ml 5M NaCl (0.2 M)

-20 ml 10% SDS (0.02%)

-100 ml 1M Tris HCl pH 8 (0.1 M)

-10 ml 0.5M EDTA (5 mM)

-830 ml ultra-pure (MilliQ) water

\*10 mg /ml Proteinase K solution: dissolve 100 mg in 5 ml MilliQ water and 5 ml 100% ethanol.

## REAGENTS

100 mg Proteinase K 3115879001 by [Sigma-aldrich](#)

### **Step 3.**

Incubate samples overnight at 55°C.

### **Step 4.**

Vortex the samples.

### **Step 5.**

Centrifuge samples for 10 min at 13000 rpm at room temperature.

### **Step 6.**

Transfer the supernatants with DNA in tubes filled with 600 µl isopropanol. Leave until all tubes are transferred.

### **Step 7.**

Mix one tube slowly at a time by inversion, DNA fibers will appear in the tube.

### **Step 8.**

Collect DNA fibers with pipette tips (put your finger on the end of the tip to avoid as much as possible that liquid enters in the tip)

### **Step 9.**

Transfer pipett tips in new labeled tubes containing 200 µl TE buffer (allow rehydration of DNA to occur during 15 min before removing the tips)

\*TE buffer (1L):

-10 ml 1 M Tris-HCl pH 7.5 (10 mM)

-2 ml 0.5 M EDTA (1 mM)

-988 ml MilliQ water

-autoclave

### **Step 10.**

Incubate 1-2 hours at 55°C and vortex to dissolve the DNA. If genotyping is postponed, store samples at -20°C, otherwise proceed.

### **Step 11.**

Prepare samples by combining PCR beads with 5µl of forward primer 5'-TTT GAA TTT GTC TTT GGA AGG CCT C-3' (10 µM work concentration; 2 µM final concentration) and 5µl of reverse primer 5'-TAT CCC CGG ATG AGT CAA CAC TGT C-3' (10 µM work concentration; 2 µM final concentration), 14µl MilliQ water and 1µl DNA.

Include negative control (no DNA).

Include DNA samples of a previous genotyped ADAMTS5-P wild-type (WT), heterozygous (HE) and homozygous (knockout, KO) animal as positive controls.

Primer sequences were provided by the group of Prof. J. Sandy (Rush University Medical Center, Chicago, IL, USA).



## REAGENTS

Illustra PuReTaq Ready-To-Go PCR Beads 27-9557-02 by [Ge Healthcare](#)

### Step 12.

Put PCR samples in 2720 Thermal Cycler (Applied Biosystems, Foster City, CA, USA)

### Step 13.

Cycling conditions are:-94°C 2 min-40 cycles of 94°C 30 sec; 63°C 30 sec; 68°C 95 sec-72°C 5 min

### Step 14.

Pour a 2% (= 2g/100 ml) agarose gel in 1x TAE buffer and immerse gel in 0.5x TAE supplemented with Midori Green Advanced DNA Stain (5µl/100 ml buffer)

\* 50x TAE (1L):

-242g Tris base

-57.1 ml glacial acetic acid

-100 ml 500 mM EDTA pH 8.0 (50 mM)

-Add MilliQ water until 800 ml

-Stir

-Bring final volume to 1L with MilliQ water

-Autoclave

\*1x TAE (1L):

-20 ml 50x TAE

-980 ml MilliQ water

Contains 40 mM Tris, 20 mM glacial acetic acid and 1 mM EDTA.

\*0.5x TAE (1L):

-10 ml 50x TAE

-990 ml MilliQ water

Contains 20 mM Tris, 10 mM glacial acetic acid and 0.5 mM EDTA.



#### REAGENTS

- ✓ Molecular Biology Grade Agarose EP-0010-05 by Contributed by users
- ✓ Midori Green Advanced DNA stain NG MG04 by Contributed by users

#### Step 15.

Prepare samples: add 5µl Orange DNA loading dye to 25µl PCR mix



#### REAGENTS


Orange DNA loading dye R0631 by [Thermo Fisher Scientific](#)

#### Step 16.

Load samples (18µl) together with 18µl TriDye 100 bp DNA ladder



#### REAGENTS

 TriDye 100 bp DNA Ladder - 125 gel lanes [N3271S](#) by [New England Biolabs](#)

#### Step 17.

Run samples for 30 min with 135V on a Mupid-One electrophoresis System Eurogentec SA

#### Step 18.

Image gel using the Biorad Molecular Imager Gel Doc XR+ Image System, take picture and digitally save the image.

#### Step 19.

Identify ADAMTS5 WT samples as samples exhibiting a single band of 650 bp, ADAMTS5 KO samples as samples showing a single band of 380 bp, while the ADAMTS5 HE samples show both bands.

## Warnings

-Dispose of Midori Green contaminated buffer and gels according to the institutional safety guidelines.