

# Colorimetric sulfo-phospho-vanillin (SPV) method for analysis of lipids in mucin Version 2

Xueying Zhong

## Abstract

**Citation:** Xueying Zhong Colorimetric sulfo-phospho-vanillin (SPV) method for analysis of lipids in mucin. **protocols.io** dx.doi.org/10.17504/protocols.io.euybexw

**Published:** 12 Apr 2016

## Protocol

### Step 1.

#### Prepare the mucin samples

Dissolve the mucins in water at a predetermined concentration;

### Step 2.

Extraction of lipids in mucin sample solution

1. For each 1 ml of sample solution, add 3.75 ml 1:2 (v/v) chloroform: methanol and vortex well;
2. Then add 1.25 ml chloroform and vortex well;
3. Add 1.25 ml dH<sub>2</sub>O and vortex well;
4. Centrifuge at 1000 RPM for 5 min to give a two-phase system (aqueous top, organic bottom);
5. Recover the bottom phase with a Pasteur pipette and store in tubes.

### Step 3.

#### Prepare the standard sample solution

1. Prepare the solvent, chloroform: methanol=2:1;
2. Mix cholesterol in solvent at predetermined concentration, for instance 5mg/ml;
3. Vary volume of the standard sample (1 ul, 2 ul, 4 ul, 8 ul, 16 ul and 24 ul) to assign different amount of cholesterol in different tubes.

### Step 4.

Evaporating solvent in both standard samples and mucin samples

### Step 5.

#### Measure background absorbance

1. Add 100 ul concentrated sulfuric acid into each tube and incubating at 90 C for 10 min ( on a dry heating bath);
2. Cooling to room temperature and measuring background absorbance at 540nm;

### Step 6.

#### Measure the absorbance after color development

1. Prepare the sulfo-phosphoric-vanillin acid agent: 0.2 mg vanillin per ml 17% phosphoric acid) for color development;
2. Add 50 ul sulfo-phosphoric-vanillin acid agent for color development;
3. Measuring absorbance at 540 nm after 5 min of color development.