



Sep 25, 2019

Q5 Polymerase PCR Cloning

NUS iGEM¹¹National University of Singapore

1

Works for me

[dx.doi.org/10.17504/protocols.io.7pqhmmw](https://doi.org/10.17504/protocols.io.7pqhmmw)

NUS iGEM

National University of Singapore



MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Q5 High-Fidelity PCR Kit - 200 rxns

E0555L

New England Biolabs

STEPS MATERIALS

NAME ▾

CATALOG # ▾

VENDOR ▾

Water refers to sterilized deionized water

Q5 Reaction Buffer Pack - 6.0 ml

B9027S

New England Biolabs

dNTP Mix, PCR Grade (200 µl)





201900

Qiagen

Q5 High-Fidelity DNA Polymerase - 100 units

M0491S

New England Biolabs

- 1 Prepare 6-reaction PCR mix by adding together  **244.2 µl** deionized water,  **66 µl** Q5 reaction buffer,  **6.6 µl** dNTP mix, and  **3.3 µl** Q5 DNA polymerase



Water refers to sterilized deionized water



Q5 Reaction Buffer Pack - 6.0 ml

by [New England Biolabs](#)

Catalog #: [B9027S](#)



dNTP Mix, PCR Grade (200 µl)

by [Qiagen](#)




Catalog #: [201900](#)












Q5 High-Fidelity DNA Polymerase - 100 units

by [New England Biolabs](#)

Catalog #: [M0491S](#)

- 2 Add  **1 µl** of template DNA and  **0.5 µl** of primer pairs into specified PCR tube
- 3 Aliquot  **48 µl** of PCR mix into each PCR tube
- 4 Load PCR tubes in thermocycler



Protocol is set to run at  **98 °C** for  **00:00:40** followed by  **00:00:30** for the desired annealing temperature before  **72 °C** for x amount of time (depending on fragment length). The cycle repeats for 34 rounds starting with  **00:00:10** at  **98 °C** followed by annealing and elongation, before incubating at  **72 °C** for  **00:05:00** and then  **4 °C** infinitely. In this PCR run, the elongation time and annealing temperature depends on specific fragment length and primer annealing temperature.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited