



Oct 20, 2019

## Golden Gate lvl 1/2

Vinca Seiler<sup>1</sup>, René Inckemann<sup>1</sup><sup>1</sup>iGEM Marburg 2019

1

Works for me

dx.doi.org/10.17504/protocols.io.8d4hs8w

iGEM Marburg  
iGEM Marburg 2019

## ABSTRACT

Golden Gate reaction protocol for lvl 1/2

## MATERIALS

NAME	CATALOG #	VENDOR
T4 DNA Ligase	M0202	New England Biolabs
Esp3I	R0734L	New England Biolabs
BsaI-HFv2	R3733S	New England Biolabs
10X NEB T4 DNA ligase buffer		New England Biolabs

## Pipetting scheme for assembly reaction

- 1 **0.5 µl** of each DNA insert. For an improved assembly efficiency, the amount of DNA inserts can optionally be normalized to equimolar concentrations (~ **20 fmol** each) or use **75 ng** of each insert (antibiotic resistance part should be diluted 1:10, **7.5 ng** - **10 ng**).
- 2 **1 µl** T4 DNA Ligase buffer (NEB)
- 3 **0.5 µl** T4 DNA Ligase (NEB)
- 4 **0.5 µl** BsaI-HF<sup>®</sup>v2(NEB) for **lvl 1** / EspI3 (NEB) for **lvl 2**
- 5 Water to **10 µl** . ( **20 µl** also possible) **0.5 µl** of each DNA insert

## Thermocycler Improved Protocol

- 6 **37 °C** **00:01:30**
- 7 **16 °C** **00:03:00**

8 Cycle step 6 and 7 15x

9  50 °C  00:05:00

10  80 °C  00:10:00

#### Alternative Thermocycler Troubleshoot/Overnight Protocol

11  37 °C  00:02:00




12  16 °C  00:05:00

13 Cycle steps 11 & 12 x 50

14  50 °C  00:10:00

15  80 °C  00:10:00

#### Transformation

16 Add  2 µl -  5 µl of each assembly reaction added to  50 µl competent cells.

17 Cells should be recovered for  01:00:00 (Amp) to  00:02:00 (Kan, Chloramphenicol).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited