

Long staining procedure of nuclei in *Euplotes crassus* using DAPI Version 2

Rachele Cesaroni

Abstract

Citation: Rachele Cesaroni Long staining procedure of nuclei in *Euplotes crassus* using DAPI. **protocols.io**

dx.doi.org/10.17504/protocols.io.hh7b39n

Published: 31 Mar 2017

Protocol

Step 1.

Pellet *Euplotes crassus* cells at 400 rcf for 3 minutes and remove as much supernatant as possible by pipetting.

NOTES

Estienne Swart 31 Mar 2017

Both algae and bacteria are autofluorescent. Better to have a completely starved *Euplotes crassus* culture.

Step 2.

Add 1 ml of 2% PFA in 1X PHEM or 4% PFA in 1X PBS to the cells and incubate them for 10 minutes at room temperature.

Step 3.

Pellet *Euplotes crassus* cells by centrifugation at 400 rcf for 3 minutes and remove as much supernatant as possible by pipetting.

Step 4.

Wash cells twice with 1X PBS (400 rcf for 3 minutes each time).

Step 5.

Add 1 ml of TBSTEM - 3% BSA and 0.5 μ l of DAPI (0.1 mg/ml) to the cells and stain for 10 minutes at room temperature.

Step 6.

Pellet *Euplotes crassus* cells by centrifugation at 400 rcf for 3 minutes.

Step 7.

Add 50 µl of Prolong medium.

Step 8.

Place an approx. 10 µl droplet of *Euplotes crassus* cells on a slide for observation by fluorescence microscopy.