Blue/White Screening of Bacterial Colonies X-Gal/IPTG Plates

Sean Seaver

Abstract

Blue/White Screening of Bacterial Colonies X-Gal/IPTG Plates

Citation: Sean Seaver Blue/White Screening of Bacterial Colonies X-Gal/IPTG Plates. protocols.io

dx.doi.org/10.17504/protocols.io.cjauid

Published: 03 Nov 2014

Guidelines

Dry X-Gal/IPTG-coated media in a laminar flow hood for approximately 30 minutes before use.

Plate Surface

- 1. Dry media plates in a laminar flow hood.
- 2. Add 40 μ l 100mM IPTG and 120 μ l X-Gal (20 mg/ml) to the surface of each plate and spread over entire surface.

Note: The edges of the plate are difficult to spread adequately and may give false positives. We advise picking colonies in the middle of the plate, if possible, for best results.

Before start

Prepare 20 mg/ml X-Gal solution in DMF (See X-Gal Stock Solution Procedure). For reduced DMF toxicity in media, you can alternatively make a 100 mg/ml X-Gal solution in DMF (this concentration is only stable at -20° C for ~ 1 week).

Prepare 100mM IPTG solution in dH2O (or dilute from 1M IPTG Stock Solution).

Materials

■ IPTG <u>GB-I2481C</u> by <u>P212121</u>

• X-Gal <u>GB-X4281C</u> by <u>P212121</u>

Protocol

Step 1.

Cool autoclaved growth media agar to 50°C.

Step 2.

Add 10 µl X-Gal Solution (20 mg/ml) per 1 mL of Media (or 2 µl X-Gal Solution (100 mg/ml) per 1 mL of Media).

O DURATION

00:02:00

Step 3.

Add 10 µl IPTG (100mM) per 1 mL of Media for a final concentration of 1mM.

REAGENTS

• IPTG <u>GB-I2481C</u> by <u>P212121</u>

O DURATION

00:03:00

Step 4.

Add screening antibiotic of choice (Ampicillin, Kanamycin, Carbenicillin, etc).

REAGENTS

Ampicillin Sodium Salt <u>TWA-A-301</u> by <u>P212121</u>

© DURATION

00:03:00

Step 5.

Pour plates and allow to cool to room temperature (usually at least 30 minutes) before use.

O DURATION

00:40:00

Step 6.

Spread transformed competent cells as desired.

Note: Blue/White Selection plates are generally stable for only 1 week if stored at 4°C in clear sleeves, but may be stored in the dark (or a dark sleeve) at 4°C for up to 1 month.

O DURATION

00:05:00