# DCX Immunohistochemistry Protocol Version 2

Nara Gyzely de Morais Magalhães , Cristovam Guerreiro Diniz and Cristovam Wanderley Picanço Diniz

### **Abstract**

Protocol for Immunohistochemistry free-floating sections with anti-doublecortin antibody for avian tissue.

Citation: Nara Gyzely de Morais Magalhães , Cristovam Guerreiro Diniz and Cristovam Wanderley Picanço Diniz DCX

Immunohistochemistry Protocol. **protocols.io** dx.doi.org/10.17504/protocols.io.h57b89n

Published: 24 May 2017

### **Materials**

- Boric acid View by P212121
- Ammonium Chloride View by P212121
- Sodium acetate by P212121
- Sodium Phosphate monobasic by Contributed by users

Triton X-100 T8787-50ML by Sigma Aldrich

Sodium phosphate dibasic <u>7558-79-4</u> by <u>Sigma Aldrich</u>

Normal Horse Serum S-2000 by Vector Laboratories

Doublecortin Antibody sc-8066 by Santa Cruz Biotechnology

Biotinylated Horse Anti-Goat IgG Antibody BA-9500 by Vector Laboratories

VECTASTAIN Elite ABC HRP Kit (Peroxidase, Standard) PK-6100 by Vector Laboratories

Ammonium nickel(II) sulfate hexahydrate 7785-20-8 by Sigma Aldrich

3,3'-Diaminobenzidine tetrahydrochloride <u>D5905</u> by <u>Sigma Aldrich</u>

#### **Protocol**

### Antigen Retrieval Method

#### Step 1.

To remove the excess paraformaldehyde fixative, wash with 0,1 M PBS the sections free-floating at agitation in room temperature for 3 min. (1/3)

**O** DURATION

00:03:00

### Antigen Retrieval Method

### Step 2.

To remove the excess paraformaldehyde fixative, wash with 0,1 M PBS the sections free-floating at agitation in room temperature for 3 min. (2/3)

© DURATION

00:03:00

### Antigen Retrieval Method

# Step 3.

To remove the excess paraformaldehyde fixative, wash with 0,1 M PBS the sections free-floating at agitation in room temperature for 3 min. (3/3)

**O DURATION** 

00:03:00

### Antigen Retrieval Method

### Step 4.

Incubate with 12% Boric acid (pH =9,0,70°C) in a water bath for 1 hour.

**O DURATION** 

01:00:00

NOTES

Nara Magalhaes 24 May 2017

When the temperature reaches 50°C, using a brush put the sections at the recipient and wait the temperature reaches 70°C and keep at this temperature for 1 hour.

#### Antigen Retrieval Method

### Step 5.

Remove the recipients from the water bath and wait until it is at room temperature.

### Antigen Retrieval Method

### Step 6.

Wash the sections with 0,1% PBS/T for 5 minutes. (1/3)

© DURATION

00:05:00

### Antigen Retrieval Method

# Step 7.

Wash the sections with 0,1% PBS/T for 5 minutes. (2/3)

**O** DURATION

00:05:00

### Antigen Retrieval Method

Step 8.

Wash the sections with 0,1% PBS/T for 5 minutes. (3/3)

**O DURATION** 

00:05:00

### Antigen Retrieval Method

Step 9.

Wash the sections with 0,1M PBS for 2 minutes (shaking). (1/3)

**O DURATION** 

00:02:00

### Antigen Retrieval Method

Step 10.

Wash the sections with 0,1M PBS for 2 minutes (shaking). (2/3)

**O DURATION** 

00:02:00

### Antigen Retrieval Method

**Step 11.** 

Wash the sections with 0,1M PBS for 2 minutes (shaking). (3/3)

**O DURATION** 

00:02:00

### **Protein Blocking Step**

Step 12.

Incubate with blocking buffer (Normal Horse Serum 10% in 0,3% PBST ) for 12 hours at gentle agitation in refrigeration 4°C.

© DURATION

12:00:00

### **Primary Antibody**

Step 13.

Remove the serum and incubate with Anti Doublecourtin antibody (Doublecortin C-18, sc-8066 Santa Cruz Biotechnology – 1:500), diluted in 0,3%PBS/T at gentle agitation, overnight at 4°C.

# **Primary Antibody**

Step 14.

Wash the sections with 0,1M PBS/T 0,1% for 2 minutes (shaking). (1/3)

**O DURATION** 

00:02:00

### **Primary Antibody**

Step 15.

Wash the sections with 0,1M PBS/T 0,1% for 2 minutes (shaking). (2/3)

**O DURATION** 

00:02:00

### **Primary Antibody**

### **Step 16.**

Wash the sections with 0,1M PBS/T 0,1% for 2 minutes (shaking). (3/3)

**O DURATION** 

00:02:00

### **Primary Antibody**

# Step 17.

Incubate with secondary antibody (Biotinylated Horse Anti-Goat IgG Antibody, BA-9500, Vector Laboratories) Diluted in 0,3% PBS/T at 1:400, for 1 hour at room temperature.

**O DURATION** 

01:00:00

### **Blocking Step**

### **Step 18.**

Incubate with 0,3% hydrogen peroxide (diluted in 0,1 M PBS) for 15 minutes with light shaking.

**O DURATION** 

00:01:00

# **Blocking Step**

### Step 19.

Wash sections with PBS/T 0,1% for 2 minutes (shaking). (1/3)

**O DURATION** 

00:02:00

### **Blocking Step**

### Step 20.

Wash sections with PBS/T 0,1% for 2 minutes (shaking). (2/3)

© DURATION

00:02:00

## **Blocking Step**

### **Step 21.**

Wash sections with PBS/T 0,1% for 2 minutes (shaking). (3/3)

**O** DURATION

00:02:00

# ABC

### Step 22.

Incubate in VECTASTAIN® ABC KIT solution for 1 hour (total) at 4°C with light shaking.

First 37,5  $\mu$ l A + 37,5 $\mu$ l B with 1,88 ml 0,3%PBS/T for 30 minutes.

© DURATION

00:30:00

#### ABC

### Step 23.

Add 13,12 ml 0,3%PBS/T and incubate for 30 more minutes at 4°C with light shaking.

© DURATION

00:02:00

### ABC

### Step 24.

Wash with 0,1% PBS /T for 5 minutes, (shaking). (1/2)

**O DURATION** 

00:05:00

#### ABC

#### Step 25.

Wash with 0,1% PBS /T for 5 minutes, (shaking). (2/2)

**O** DURATION

00:05:00

### DAB Visualization - GDN preparation

#### Step 26.

Firstly prepare the Solution A by mix 0,006g of Diaminobenzidine (DAB) with 5 ml of distilled water.

### DAB Visualization - GDN preparation

### **Step 27.**

Secondly, prepare the Solution B by mixing 0,250g of Nickel ammonium sulfate with 5 ml de Acetate Buffer pH 6.0.

# DAB Visualization - GDN preparation

### **Step 28.**

Thirdly mix Solution A and B adding ammonium chloride (0,004g) with 0.020g  $\alpha$ -D-Glucose.

### DAB Visualization - GDN preparation

#### Step 29.

Leave the section in this mix of A and B for 5 minutes.

### © DURATION

00:05:00

#### **P** NOTES

# Nara Magalhaes 24 May 2017

After each step remove the remaining solution from the previous step with a pipette, with careful not to damage or lose sections.

Use only sterilized material to minimize the risk of contamination between different antibodies or solutions.

At the end add chlorine to inactivate GND + DAB residual solution before disposal on an appropriate container.

### **DAB Visualization**

#### Step 30.

Incubate sections with solution GND and wait for 3 minutes.

**O DURATION** 

00:03:00

### **DAB Visualization**

### **Step 31.**

Add 0,007g of Glucose-oxidase for each 3ml of GND solution for revelation. Stop revelation when the goal contrast is achieved (use a low gain microscope).

# **DAB Visualization**

#### Step 32.

Remove the GDN + Glucose oxidase and wash the sections using 0,1M PBS for 2 minutes with light shaking. (1/3)

© DURATION

00:02:00

### **DAB Visualization**

# **Step 33.**

Remove the GDN + Glucose oxidase and wash the sections using 0,1M PBS for 2 minutes with light shaking. (2/3)

© DURATION

00:02:00

### **DAB Visualization**

## **Step 34.**

Remove the GDN + Glucose oxidase and wash the sections using 0,1M PBS for 2 minutes with light

shaking. (3/3)

© DURATION

00:02:00

# **DAB Visualization**

# Step 35.

Mount the sections in appropriate gelatinized microscope slides and dry at room temperature for 12 hours or more depending on the mounting medium of choice.

© DURATION

01:00:00

### **DAB Visualization**

# **Step 36.**

Dihydrate and add the coverslips.