

# Immunocytochemistry Staining Protocol Version 3

BioLegend, Inc., Kelsey Knight

## Abstract

**Citation:** BioLegend, Inc., Kelsey Knight Immunocytochemistry Staining Protocol. **protocols.io**

dx.doi.org/10.17504/protocols.io.tdmei46

**Published:** 10 Sep 2018

## Before start

### Reagent List:

- Chamber slides, cover slips, or 12-well plates
- Phosphate-buffered saline (PBS)
- Fixation solution: 1% Paraformaldehyde, in PBS
- Permeabilization solution: 0.5% Triton X-100 in PBS
- Blocking buffer: 5% FBS in PBS

## Protocol

### Sterilization (for 12-well plates with coverslips)

#### Step 1.

Transfer a single cover slip into a 12-well plate. Then add 1mL of 70% Ethanol into a well for 20 minutes at room temperature.

 **DURATION**

00:20:00 :

### Sterilization (for 12-well plates with coverslips)

#### Step 2.

Wash quickly three times with PBS.

### Poly-Lysine Coating for 12-Well Plates (optional; for loosely attached cells)

#### Step 3.

Add 1 mL of 0.1 mg/mL Poly-D-lysine solution into a well for 15 minutes at room temperature.

 **DURATION**

00:15:00 :

### Poly-Lysine Coating for 12-Well Plates (optional; for loosely attached cells)

#### Step 4.

Wash quickly three times with PBS and let dry before plating cells.

#### Sample Preparation

##### **Step 5.**

Grow cultured cells on cover slips or in wells overnight at 37°C. At the time of fixation, cells should be 70-80% confluent in single layer.

#### Sample Preparation

##### **Step 6.**

Rinse cells briefly in PBS.

#### Sample Preparation

##### **Step 7.**

Fix cells by incubation with freshly made 1% Paraformaldehyde in PBS for 10 minutes at room temperature.

 DURATION

00:10:00 :

#### Sample Preparation

##### **Step 8.**

Rinse three times quickly in PBS.

#### Sample Preparation

##### **Step 9.**

For intracellular staining, add permeabilization solution and incubate at room temperature for 10 minutes. Then wash quickly three times in PBS.

 DURATION

00:10:00 :

#### Sample Blocking

##### **Step 10.**

Block samples in 1 mL of blocking buffer at room temperature for 30 minutes.

 DURATION

00:30:00 :

#### Sample Staining

##### **Step 11.**

Dilute the primary antibody to the recommended concentration/dilution in blocking buffer.

#### Sample Staining

##### **Step 12.**

For 8-well chamber slides, add 200 µL per well. For 12-well plates, add 500 µL per well. Incubate two

to three hours at room temperature or overnight at 4°C. If using conjugated antibodies, perform this step in the dark.

### Sample Staining

#### Step 13.

For surface staining, rinse 3 times quickly in PBS. For intracellular staining, quickly wash once followed by incubation with wash buffer for 5-10 minutes. Then quickly wash additional two times.

**Note:** If using primary antibodies directly conjugated to fluorochromes, then skip to step 17.

 DURATION

00:10:00 :

### Sample Staining

#### Step 14.

Prepare fluorochrome-conjugated secondary antibody in blocking buffer according to the manufacturer's specification data sheet, and add 200 µl per well to the 8-well chamber slides. For 12-well plates, add 500 µL per well.

### Sample Staining

#### Step 15.

Incubate the samples for one hour, at room temperature, in the dark.

 DURATION

01:00:00 :

### Sample Staining

#### Step 16.

For surface staining, rinse three times quickly in PBS. For intracellular staining, quickly wash once followed by incubation with wash buffer for 5-10 minutes, then quickly wash additional two times.

 DURATION

00:10:00 :

### Sample Staining

#### Step 17.

Optional: To stain F-actin, prepare a working solution of Flash Phalloidin™ by diluting it 1:20-1:100 in PBS. Add 200 µL per well for an 8-well plate or 500 µL per well for a 12-well plate. Stain for 20 minutes at room temperature in the dark.

 DURATION

00:20:00 :

### Sample Staining

#### Step 18.

Apply anti-fade mounting medium to the cover slip.

**Step 19.**

Seal slides with nail polish.