



Aug 07,
2019

Quantitative PCR analysis to assess gene expression changes in hyperglycemic larval zebrafish



Sandra Rieger¹

¹University of Miami

1

Works for me

dx.doi.org/10.17504/protocols.io.3bagiie

Diabetic Complications Consortium

Tech. support email: rmcindoe@augusta.edu



Lili Liang



ABSTRACT

Summary:

To assess hyperglycemia in larval zebrafish, we analyzed the expression of genes implicated in glucose metabolism, such as insulin, insulin receptor, glucagon, and phosphoenolpyruvate carboxylase (pepck), using quantitative PCR (qPCR). We analyzed gene expression following treatment of Tg(*ins:NTRmCherry*) transgenic larvae either with 0.5 % DMSO (controls) or 10 mM metronidazole (see Protocol 2) for β -cell ablation at 3 and 8 days post fertilization (dpf).

Diabetic Complication:



Neuropathy

EXTERNAL LINK

<https://www.diacomp.org/shared/document.aspx?id=222&docType=Protocol>

MATERIALS

NAME	CATALOG #	VENDOR
RNeasy Mini Kit	74104	Qiagen
oligo (dT) or random hexamer primers		Integrated DNA Technology
SuperScript [®] III First-Strand Synthesis System for RT-PCR	18080-051	Life Technologies
SYBR Green Brilliant III Ultra-Fast SYBR [®] Green QPCR Master Mix	600882	Agilent Technologies
ROX reference dye (included in Brilliant III master mix)	600882	Agilent Technologies
Mx3000P QPCR System	Mx3000P	Stratagene

MATERIALS TEXT

Primers:

Gene name	Gene symbol	5'-3' forward primer	5'-3' reverse primer
<i>insulin</i>	<i>ins</i>	ccc ttt atc tgg tct gtg gc	ttc ctt atc agc tcg gca tg
phosphoenolpyruvate carboxylase	<i>pepck</i>	gag aac agc acc atc ctc ag	tca ccg ttt tac tct cca cac
<i>glucagon a</i>	<i>gcga</i>	aag act tcg ttc agt ggc tc	tgg gtt gtc cgg att tta gc
<i>insulin receptor a</i>	<i>insra</i>	ata aga atg atc ggg agt gtg g	tgc ttg cat gag gat gga c
<i>Elongation factor alpha (control gene)</i>	<i>efla</i>	tct aca aat gcg gtg gaa tcg	gag caa tgt caa tgg tga tac c

Note:

QIAGEN ([RRID:SCR_008539](#))

- 1 Preparation of mRNA from larval zebrafish. Pool 10 zebrafish larvae and isolate RNA according to the RNeasy Mini Kit (QIAGEN) manual. Elute final mRNA in 30 µl of RNase-free water.
- 2 Preparation of cDNA from mRNA. Prepare cDNA from isolated mRNA according to the SuperScript®III First-Strand Synthesis System (Life Technologies) using either oligo(dT) or random hexamer primers.
- 3 Preparation of qPCR mix (according to the Agilent manual for preparation of the Brilliant III Ultra- Fast SYBR® Green QPCR Master Mix):

Dilute the reference dye 1:500 using nuclease-free PCR-grade water. Prepare the experimental reactions by combining the components of the reagent mixture in the order listed in the table below. Prepare a single reagent mixture for replicate reactions (plus at least one reaction volume excess) using multiples of each component.

Reagent Mixture

Nuclease-free PCR-grade water to bring final volume to 20 µl (including cDNA)
 10 µl of 2× SYBR Green QPCR Master Mix
 x µl of upstream primer at optimized concentration (200–500 nM)
 x µl of downstream primer at optimized concentration (200–500 nM)
 0.3 µl of diluted ROX reference dye

Gently mix the reagent mixture without creating bubbles, then distribute the mixture to the experimental reaction tubes. Add 0.5-50 ng (x µl) of cDNA to each reaction to bring the final reaction volume to 20 µl.

We amplified each transcript using the following PCR conditions:

- 1) 95°C 3 minutes
- 2) **40 cycles**
 95°C 15 seconds
 60°C 20 seconds,
- 3) **Melting curve**
 95°C 1 minute
 60°C 30 seconds
 95°C 30 seconds

Potential Pitfalls:

- 4 Inconsistent qPCR results: Always use master mixes and the suggested reference dye (ROX reference dye, Life Technologies, USA). Include pre-mix of triplicate reactions with sample and primers. Alternatively, utilize the Taqman system (Life Technologies).



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited