

# WB transfer buffer (1L of 10X)

Bao Thai

## Abstract

**Citation:** Bao Thai WB transfer buffer (1L of 10X). **protocols.io**

[dx.doi.org/10.17504/protocols.io.mpac5ie](https://dx.doi.org/10.17504/protocols.io.mpac5ie)

**Published:** 16 Jan 2018

## Protocol

### Step 1.

10X WBTB (1L)  
250 mM Tris pH 8.0 (30.2 g)  
1.92 M Glycine (144 g)

1X Transfer Buffer (TB; 1L)  
(make day of)  
100 ml 10X WBTB  
200 ml MeOH  
700 ml cold water