

# Plaque Assay Protocol for Pseudoalteromonas

Natalie Solonenko

## Abstract

**Citation:** Natalie Solonenko Plaque Assay Protocol for Pseudoalteromonas. **protocols.io**

[dx.doi.org/10.17504/protocols.io.de33gm](https://dx.doi.org/10.17504/protocols.io.de33gm)

**Published:** 02 Feb 2016

## Before start

### What you need before you start:

1. 20% nutrient Zobell plates
2. Top agar – 3.5 ml per plate
  - a. 100% nutrient Zobell
  - b. 6g agar/liter
3. Your host growing somewhere in exponential phase – 0.4 ml per plate
4. Phages – this could be:
  - a. A lysate
  - b. A plaque picked into buffer
  - c. An environmental sample
  - d. A sample from an experiment (e.g., a one-step)
5. A 40°C water bath

## Protocol

### Step 1.

First, autoclave the top agar to melt it

### Step 2.

While the autoclave is running:

- a. Let your plates warm to room temperature and label them
- b. Make whatever dilutions of your phages you plan to plate

### Step 3.

Aliquot the agar into tubes and place in 40°C water bath

### NOTES

**VERVE Team** 15 Jul 2015

Before you use the agar, make sure the tube feels the same temperature as the water bath. If the agar is still too hot, you'll kill your host. You'll need to wait about 5–15 minutes for it to cool down enough. If the agar is too cool, it might start getting chunky. If your agar is starting to solidify, do NOT use it. Make up a new tube if your big bottle is still hot or re-autoclave it. If it will take you more than 30 minutes to do all your infections, you're doing too many at a time – your agar will start to solidify.

**VERVE Team** 15 Jul 2015

If you're doing more than one plate with the same host, you can aliquot enough agar for however many plates you're comfortable doing at once (up to 12) Be sure to leave room in your tubes to add host.

**VERVE Team** 15 Jul 2015

Use 35 ml per plate.

#### **Step 4.**

Pipet 100 µl of your phage sample onto the plate

#### **NOTES**

**VERVE Team** 15 Jul 2015

If you made up agar aliquots for more than one plate, add phage to each plate now.

#### **Step 5.**

Remove the agar tube from the water bath and add the host to it

#### **NOTES**

**VERVE Team** 15 Jul 2015

Pipet up and down or gently invert to mix. Try not to introduce a lot of bubbles (i.e., do not shake or vortex the tube).

**VERVE Team** 15 Jul 2015

Use 0.4 ml of host per plate.

#### **Step 6.**

Add 3.5 ml of the agar/host mixture

#### **NOTES**

**VERVE Team** 15 Jul 2015

If you want your plates in stacks, stack them as you go. They cannot be moved after plating for at least an hour.

**VERVE Team** 15 Jul 2015

If there is only one plate worth of agar and host in the tube, you can pour it onto the plate. Otherwise, pipet 3.5 ml onto each plate.

#### **Step 7.**

Swirl each plate to spread the agar as you go

#### **NOTES**

**VERVE Team** 24 Aug 2015

Make sure it covers the whole plate.

#### **Step 8.**

Leave the plates on the bench overnight

#### **DURATION**

18:00:00

#### **Step 9.**

Count or pick plaques the next day

#### **NOTES**

**VERVE Team** 15 Jul 2015

If you're calculating PFU, you will probably need to leave them out and count them for one more day.