

# Staining protocol for unicellular protists: LysoTracker and Dapi

Maria Rubio-Brotons

## Abstract

**Citation:** Maria Rubio-Brotons Staining protocol for unicellular protists: LysoTracker and Dapi. **protocols.io**  
dx.doi.org/10.17504/protocols.io.geybtfw

**Published:** 16 Nov 2016

## Protocol

### 2% PFA Fixation

#### Step 1.

Add PFA drop by drop and mixing carefully into 1mL cells for a final concentration of 2%PFA

Incubate cells for 4 min at room temperature

#### Step 2.

Centrifuge at 1000g for 5 minutes

#### Step 3.

Wash cells with 1XPBS twice

#### Step 4.

Resuspend cells with 1mL 1xPBS

#### Step 5.

Add 1/1000 lysoTracker dye

#### Step 6.

1 uL of lysoTracker for a 1mL final Volume of cells

Room temperature, protect from light

After 20 minutes, add 1/1000 DAPI dye

#### Step 7.

Room temperature, protect from light

After 5 minutes, wash with 1xPBS twice

#### Step 8.

Centrifuge at 1000g for 5 minutes and resuspend with 1xPBS

Total incubation times:

Lysotracker: 25min

DAPI: 5min

Mount your slide and ready for observation!

**Step 9.**