



Zymo Duet DNA/RNA MiniPrep Plus Extractions

Amy Lyden¹

¹CZ Biohub

[dx.doi.org/10.17504/protocols.io.umteu6n](https://doi.org/10.17504/protocols.io.umteu6n)

Amy Lyden

ABSTRACT

Attached are detailed instructions for following the kit guidelines for the Zymo Duet DNA/RNA MiniPrep Plus.

[_D7003_ZR-Duet_DNA-RNA_MiniPrep_Plus_ver_1.0.1.pdf](#)

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

GUIDELINES

Follow kit guidelines for best practices.

[_D7003_ZR-Duet_DNA-RNA_MiniPrep_Plus_ver_1.0.1.pdf](#)

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
ZR-Duet™ DNA/RNA MiniPrep Plus111212	D7003	Zymo Research

MATERIALS TEXT

Reagent in Kit	Quantity	Storage Temperature	Notes on Preparation
DNA/RNA Lysis Buffer	50 ml	Room Temp.	
DNA/RNA Prep Buffer	50 ml	Room Temp	
DNA/RNA Wash Buffer (concentrate)	24 ml	Room Temp	Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml DNA/RNA Wash Buffer concentrate before use
DNase/RNase-Free Water	10 ml	Room Temp	
DNase I (lyophilized)	1	Room Temp (lyophilized), -20 C (reconstituted)	Prior to use, reconstitute the lyophilized DNase I with 275 µl DNase/RNase-Free Water. Mix by inversion. Store aliquots at -20°C.
DNA Digestion Buffer	4 ml	Room Temp	
DNA/RNA Shield™ (2x concentrate)	25 ml	Room Temp	Dilute in DNase/RNase-free Water to obtain 1X for storage
PK Digestion Buffer	5 ml	Room Temp	
Proteinase K	20 mg	Room Temp (lyophilized), -20 C (reconstituted)	Prior to use, reconstitute the lyophilized Proteinase K with 1040 µl Proteinase K Storage Buffer. Vortex to dissolve. Store at -20°C
Spin-Away™ Filters	50	Room Temp	
Zymo-Spin™ IIIICG Columns	50	Room Temp	
Collection Tubes	3x50	Room Temp	
Other Reagents and Equipment			
Bead bashing tubes (ZR BashingBead™ Lysis Tubes or Omni 2.8mm ceramic bead bashing tubes or other)			
Qiagen Tissulyser II			
Filtered Pipette tips			
Ethanol			
Eppendorfs (1.5 mL)			
Microcentrifuge			
Eppendorf rack			





SAFETY WARNINGS

Operate in an appropriate BSL level space for your sample type and suspected infectious agents.





BEFORE STARTING

Ensure your samples have been stored at -80C, preferably in RNA/DNA shield. Avoid freeze thaws to ensure sample integrity. Identify your sample type and follow the 'sample preparation' section for that type.



Sample Preparation for cultured isolates

1. Add  **320 µl** 1X DNA/RNA shield (2x concentrate diluted in nuclease-free water) to bead bashing tubes
2. Select colonies and add to DNA/RNA shield in bead bashing tube.
3. Tissue lyse on Qiagen Tissue lyser II  **00:02:00 minutes** at a frequency of 1/30 s. Alternatively, vortex at the highest frequency for 2 minutes.
4. Spin  **00:01:00** at 14000RCF.
5. Add  **320 µl** of DNA/RNA Lysis buffer. Proceed with Nucleic Acid Purification.






Sample Preparation for liquids (ie, plasma, tracheal aspirate, CSF)

2. 1. Add 1 part sample to 3 parts DNA/RNA lysis buffer (ie,  **100 µl** sample to  **300 µl** lysis buffer) in a bead bashing tube
2. Tissue lyse on Qiagen Tissue lyser II  **00:01:00 minutes** at a frequency of 1/30 s. Alternatively, vortex at the highest frequency for 1 minutes.
3. Spin  **00:01:00** at 14000RCF. Proceed with Nucleic Acid Purification.




Sample Preparation for samples in DNA/RNA shield (1x)




3. 1. If not homogenized, tissue lyse on Qiagen Tissue lyser II in bead bashing tube for desired amount of time at a frequency of 1/30 s. Alternatively, vortex at the highest frequency.
2. Add one part DNA/RNA Lysis buffer (ie  **200 µl** of sample in DNA/RNA shield to  **200 µl** of lysis buffer). Proceed with Nucleic Acid Purification.

Sample Preparation for whole blood

4. 1. Add  **200 µl** of 2X DNA/RNA shield to each  **200 µl** of b2 of blood sample and mix thoroughly.
2. For every  **400 µl** of reagent/blood mixture, add  **8 µl** of reconstituted Proteinase K and mix thoroughly.
3. Incubate at at room temperature for  **00:30:00 minutes**
4. Add an equal volume of isopropanol to the reagent/blood mixture and mix by vortexing. Proceed with Nucleic Acid Purification.

Nucleic Acid Purification

5. Label a Spin-Away Filter (yellow) and the microcentrifuge tube for each sample.
6. Transfer up to  **700 µl** of the sample into a Spin-Away™ Filter (yellow) in a Collection Tube and centrifuge  **00:00:30 seconds** at 16000 RCF. If you have >  **700 µl** , you can reload the column and centrifuge again. Save the flow-through if you wish to proceed with RNA Purification. If not, discard flow through, and skip RNA Purification steps.

(**For whole blood ONLY**, after doing the above, discard flow-through. Add  **200 µl** of DNA/RNA Lysis buffer directly to the filter matrix. Let stand for  **00:05:00 minutes** . Then centrifuge  **00:00:30 seconds** at 16000 RCF. **SAVE** this flow through as your RNA is in it)
7. Proceed using the yellow column with DNA Purification section, and using the flow-through with the RNA purification section. After that, you will continue with the Nucleic Acid Purification for both DNA and RNA in parallel.

DNA Purification

8. DNA is bound to the yellow column. Transfer the Spin-Away™ Filter (yellow) into a new Collection Tube. Proceed with Nucleic Acid Purification.

RNA Purification

- 9 Transfer flow-through into Eppendorf if necessary. Add an equal volume of ethanol (95-100%) to the flow-through and mix well.
- 10 Transfer up to **700 µl** of the sample with ethanol into a Zymo-Spin™ III CG Column (green) in a Collection Tube and centrifuge **00:00:30 seconds** at 16000 RCF. Discard the flowthrough. Repeat by reloading the column if you have more sample remaining.
- 11 Wash the column with **400 µl** DNA/RNA Wash Buffer and centrifuge **00:00:30 seconds** at 16000 RCF. Discard the flowthrough.
- 12 Add 80 µl DNase I Reaction Mix (**5 µl** reconstituted DNase + **75 µl** DNA Digestion Buffer) directly to the column matrix. (can make a master mix of the same proportions for multiple samples)
- 13 Incubate the column at room temperature (**20 °C** - **30 °C**) for **00:15:00 minutes** . Proceed with Nucleic Acid Purification.

Nucleic Acid Purification

- 14 Add **400 µl** DNA/RNA Prep Buffer to your (yellow/green) column and centrifuge **00:00:30 seconds** at 16000RCF. Discard the flow-through.
- 15 Add **700 µl** DNA/RNA Wash Buffer and centrifuge **00:00:30 seconds** at 16000RCF. Discard the flow-through.
- 16 Add **400 µl** DNA/RNA Wash Buffer and centrifuge the column for **00:02:00 minutes** at 16000RCF to ensure complete removal of the wash buffer. Carefully transfer the column into clean microcentrifuge/Eppendorf tubes.
- 17 Add **50 µl** DNase/RNase-Free Water directly to the column matrix, let stand for **00:05:00 minutes** , and then centrifuge **00:00:30 seconds** at 16000RCF to elute DNA and RNA from the respective column.
- 18 Take all **50 µl** of eluted DNA/RNA in microcentrifuge tube and reload onto column. Let stand for **00:05:00 minutes** , and then centrifuge **00:01:00 minute** at 16000RCF to elute more DNA and RNA from the respective column.
- 19 Aliquot RNA into single use aliquot tubes and freeze at -80. Freeze DNA at -20 or -80.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited