

# In vitro digestion of DNA with Cas9 Nuclease, *S. pyogenes* (M0386) Version 3

Breton Hornblower

## Abstract

Cas9 Nuclease, *S. pyogenes* (Cas9), is a double-stranded DNA endonuclease that is guided to its target by sequence complementarity of a small RNA loaded into the protein. This protocol describes how to digest double-stranded DNA in vitro using Cas9 and a single guide RNA (sgRNA).

**Citation:** Breton Hornblower In vitro digestion of DNA with Cas9 Nuclease, *S. pyogenes* (M0386). protocols.io dx.doi.org/10.17504/protocols.io.idhca36

**Published:** 08 Jun 2017

## Guidelines

### OVERVIEW:

Cas9 Nuclease, *S. pyogenes*, (Cas9) is a double-stranded DNA endonuclease that is guided to its target by sequence complementarity of a small RNA loaded into the protein. This protocol describes how to digest double-stranded DNA in vitro using Cas9 and a single guide RNA (sgRNA).

### REQUIRED MATERIALS:

- Cas9 Nuclease, *S. pyogenes* (NEB [#M0386](#))
- 10X Cas9 Nuclease Reaction Buffer
- Nuclease-free water
- sgRNA containing the targeting sequence in the region of interest
- sgRNAs can be generated by in vitro transcription using the HiScribe T7 Quick High-Yield RNA synthesis Kit (NEB [#E2050](#)) using linearized plasmid, PCR products, or oligonucleotides as templates
- sgRNAs must contain sequence complementary to the target DNA (1,2) For information on design of sgRNA transcription templates please visit [Addgene](#)
- DNA substrate containing the target sequence
- The substrate DNA can be circular or linearized plasmid, PCR products, or synthesized oligonucleotides
- Proteinase K, Molecular Biology Grade (NEB [#P8107S](#))

### OPTIONAL MATERIALS:

Apparatus and reagents for DNA fragment analysis

- E. g. Agarose gel electrophoresis apparatus
- DNA Loading Dye (e.g. Gel Loading Dye, Purple (6X) NEB [#B7024S](#))
- E.g. Agilent Bioanalyzer or similar

## BEFORE YOU START:

We strongly recommend wearing gloves and using nuclease-free tubes and reagents to avoid RNase contamination. Further recommendations for avoiding ribonuclease contamination can be found [here](#). Reactions are typically 30 µl but can be scaled up as needed. Reactions should be assembled in nuclease-free microfuge tubes or PCR strip tubes.

It is essential to keep the molar ratio of Cas9 and sgRNA per target site at 10:10:1 or higher to obtain the best cleavage efficiency. A calculator can be found [here](#).

Prepare 300nM sgRNA by diluting the stock with nuclease-free water on ice.

Prepare 30nM substrate DNA with a single target sequence by diluting the stock with nuclease-free water on ice.

## PROCEDURE:

1. Assemble the reaction at room temperature in the following order:

Component	Volume (for 30 µl reaction)
Nuclease-free water	20 µl
10X Cas9 Nuclease Reaction Buffer	3 µl
300nM sgRNA	3 µl (30nM final)
1 µM Cas9 Nuclease, <i>S.pyogenes</i> ( <a href="#">M0386S</a> )	1 µl (~30nM final)
Reaction volume	27 µl

### Pre-incubate for 10 minutes at 25°C

30nM substrate DNA	3 µl (3nM final)
Total reaction volume	30 µl

\*The substrate DNA and sgRNA, and nuclease-free water are not included.

2. Mix thoroughly and pulse-spin in a microfuge.
3. Incubate at 37°C for 1 hour.
4. Proceed with fragment analysis.

## REFERENCES:

1. Jinek et al. (2012) Science 337 (6096) 816-821.
2. Larson et al. (2013) Nature Protocol 8 (2180-2196).
3. Mali et al. (2013) Science 339 (6121): 823-826.

## Before start

We strongly recommend wearing gloves and using nuclease-free tubes and reagents to avoid RNase contamination. Further recommendations for avoiding ribonuclease contamination can be found [here](#).

## Protocol

### Step 1.

Prepare 300nM sgRNA by diluting the stock with nuclease-free water on ice.

## Step 2.

Prepare 30nM substrate DNA with a single target sequence by diluting the stock with nuclease-free water on ice.

## Step 3.

Assemble the reaction at room temperature in the following order (total volume **30 µl**)

### **PROTOCOL**

#### **. Cas9 M0386 Mixture**

CONTACT: [New England Biolabs](#)

### **NOTES**

**New England Biolabs** 10 May 2017

Reactions are typically 30 µl but can be scaled up as needed. Reactions should be assembled in nuclease-free microfuge tubes or PCR strip tubes.

**New England Biolabs** 10 May 2017

It is essential to keep the molar ratio of Cas9 and sgRNA per target site at 10:10:1 or higher to obtain the best cleavage efficiency. A calculator can be found [here](#).

**New England Biolabs** 10 May 2017

The substrate DNA and sgRNA, and nuclease-free water are not included.

#### **Step 3.1.**

Nuclease-free water **20 µl**

#### **Step 3.2.**

10X Cas9 Nuclease Reaction Buffer **3 µl**

#### **Step 3.3.**

300nM sgRNA **3 µl** (30nM final)

#### **Step 3.4.**

1 µM Cas9 Nuclease, *S.pyogenes* (30nM final)

### **AMOUNT**

1 µl Additional info:

### **REAGENTS**

 Cas9 Nuclease, *S. pyogenes* - 70 pmol [M0386S](#) by [New England Biolabs](#)

## Step 4.

Pre-incubate for 10 minutes at 25 °C

### **DURATION**

00:10:00

## Step 5.

Add 30nM substrate DNA

## Step 6.

Mix thoroughly.

## Step 7.

Incubate at 37°C for 15 minutes.

### **DURATION**

01:00:00

**Step 8.**

Add 1  $\mu$ l of Proteinase K to each sample. Mix thoroughly and pulse-spin in a microfuge.

**Step 9.**

Incubate at room temperature for 10 minutes.

**Step 10.**

Proceed with fragment analysis.