

Radioactive Northern Blot

Anna Behle

Abstract

Protocol for Northern Blot using radioactive RNA probes.

Citation: Anna Behle Radioactive Northern Blot. **protocols.io**

dx.doi.org/10.17504/protocols.io.gs2bwge

Published: 23 Mar 2017

Before start

Generate a PCR product with a T7 promoter for generating riboprobes. Make sure the riboprobe is antisense to the RNA of interest.

Protocol

Step 1.

Separate RNA by PAGE:

✓ PROTOCOL

. [Denaturing RNA Urea-PAGE](#)

CONTACT: [Anna Behle](#)

Gel preparation

Step 1.1.

Before preparing the gel, clean all components with 70 % ethanol and RNase away.

Recipe for one small Hoefer gel (10 mL):

PAA percentage	5 %	6 %	8 %	10 %	12 %	15 %
40 % PAA (19:1)	1.25 mL	1.5 mL	2 mL	2.5 mL	3 mL	3.74 mL
Urea	5 g	5 g	5 g	5 g	5 g	5 g
10x TBE	1 mL	1 mL	1 mL	1 mL	1 mL	1 mL

10 % APS	80 µL	80 µL	80 µL	80 µL	80 µL	80 µL
TEMED	10 µL	10 µL	10 µL	10 µL	10 µL	10 µL
ddH₂O	ad 10 mL	ad 10 mL	ad 10 mL	ad 10 mL	ad 10 mL	ad 10 mL

Gel preparation

Step 1.2.

Dissolve 5 g urea in PAA, 10x TBE and H₂O. This can be done at RT or at 40°C. Agitate/shake occasionally.

SAFETY INFORMATION

Wear vinyl gloves when handling PAA (neurotoxic!)

Gel preparation

Step 1.3.

Let the liquid cool down before pouring gel.

Assemble gel caster. Take a 1 mL aliquot of gel mixture, add 20 µL APS and 1 µL TEMED and quickly pour to prevent gel from running out.

Gel preparation

Step 1.4.

Add 80 µL APS and 12 µL TEMED to the rest of the gel and quickly pour. Insert comb immediately and clasp tight to prevent leakage. Let polymerize for an hour or over night.

Gels can be stored for a week by wrapping with wet paper towels and Saran wrap.

RNA preparation

Step 1.5.

For Northern Blot Analysis, at least 2-3 µg of RNA should be loaded per well.

Adjust concentration and volume of all RNA samples by adding RNase-free ddH₂O. Add desired amount of 2x or 5x RNA loading dye.

Denature RNA at 95°C for 5 min. Snap cool on ice.

DURATION

00:05:00

SAFETY INFORMATION

RNA loading dye contains formamide. Wear goggles/lab coat/ gloves!

RNA loading

Step 1.6.

Add running buffer (1x TBE). Remove comb. Before loading samples, wash all wells and each well individually just before loading with 1xTBE and a syringe or a pipette tip. Carefully pipet samples into wells.

Step 1.7.

Close the lid and plug electrodes into power supply.

Separate RNA at 20 mA/gel, 1 hour.

Visualization

Step 1.8.

Pour out buffer, then disassemble gel.

Incubate gel in 0.5x TBE + GelRed for 10 minutes, then visualize.



00:10:00

Buffers needed

Step 2.

- Blotting Buffer: 0.5x TBE (50 mM Tris-HCl, 50 mM boric acid, 1 mM EDTA)
- 5x Ribomax buffer: 400 mM Hepes-KOH pH 7.5, 60 mM MgCl₂, 10 mM spermidine, 200 mM DTT
- 20x Roti-Stock SSC buffer (Roth): 3 M NaCl, 300 mM sodium citrate, pH=7.0
- Hybridization buffer: 50% (v/v) formamide, deionized, 7 % (w/v) SDS, 120 mM Na-Phosphate buffer, pH=7.2 (Na₂HPO₄, NaH₂PO₄, 250 mM NaCl)
- Washing buffer I: 2x SSC, 1 % (w/v) SDS
- Washing buffer II: 1x SSC, 0.5 % (w/v) SDS
- Washing buffer III: 0.5x SSC, 0.1 % (w/v) SDS

Blot

Step 3.

Measure the size of the gel.

Cut a piece of nylon membrane, as well as 4-6 Whatman papers to the corresponding size.

Incubate in 0.5xTBE.

Blot

Step 4.

On the BioRad TransBlot semi-dry blotting system, assemble blotting sandwich:

(-)Cathode

3 layers of Whatman paper

Nylon membrane

Gel

3 layers of Whatman paper

(+)Anode

Make sure to remove any air bubbles.

Blot

Step 5.

Blot at 2 mA/cm² for 45 minutes.

 DURATION

00:45:00

Crosslinking

Step 6.

Crosslink transferred RNA with the membrane using UV-light. This can be done while documenting transfer.

In vitro transcription

Step 7.

 PROTOCOL

. [Radioactive in vitro transcription](#)

CONTACT: [Anna Behle](#)

Buffers and reagents required

Step 7.1.

- 5x Ribomax buffer: 400 mM Hepes-KOH, pH 7.5; 60 mM MgCl₂; 10 mM spermidine; 200 mM DTT

Pipetting scheme

Step 7.2.

PCR Product containing the T7-Promoter (100 ng - 1 µg)

Ribomax buffer 6 µL

Pyrophosphatase 1 µL

RNase inhibitor 0.5 µL

DTT (1 M) 2 µL

ATP 10 mM 1.8 µL

CTP 10 mM 1.8 µL

GTP 10 mM 1.8 µL

α -[³² P]-UTP	3 μ L
T7 RNA-Polymerase	1 μ L
H ₂ O <i>ad</i>	30 μ L

Incubate reaction for at least 2 hours at 37°C.

 DURATION

02:00:00

Clean-up

Step 7.3.

Remove excessive radioactive nucleotides by cleaning the probe using a G-50 column (GE-Healthcare).

Before start:

Step 7.4.

Generate a PCR product using primers containing the T7-promoter. Gel-extract or clean up using column purification.

Step 8.

While the probe is transcribed:

Prehybridize membrane at 42 °C for 1 hour in hybridization buffer.

 DURATION

01:00:00

Step 9.

Add radioactive probe directly to the hybridization buffer. Make sure not to pipet onto membrane.

Hybridize for 6 hours or over night. Use calculated annealing temperature or 62°C.

Screen prepara

Step 10.

Prepare phosphor screen by illuminating for at least 30 minutes. This erases previously obtained signals.

Stringency washes:

Step 11.

Wash membrane:

- 10 min in Washing buffer I (42 °C)
- 5 min in Washing buffer II (42 °C)
- A few seconds in Washing buffer III (RT)

Detection

Step 12.

Wrap membrane in transparent plastic foil to prevent leakage of radioactive material.

Place a phosphor plate/screen on top.

Incubate depending on signal strength.

5S rRNA: 30 minutes

Other RNA: An hour to over night

Step 13.

Detect signals using a phosphoimager (Fujifilm).

Warnings

Always wear a lab coat, goggles and gloves when working with radioactive materials.

Limit exposure to radioactive materials.

Measure all materials with a Geiger counter.