



Metabolomics

 PLOS Genetics

Wen Aw¹

¹z3314717@unsw.edu.au

[dx.doi.org/10.17504/protocols.io.rtd6ie](https://doi.org/10.17504/protocols.io.rtd6ie)

Cage Studies

 Wen Aw 

EXTERNAL LINK

<https://doi.org/10.1371/journal.pgen.1007735>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Aw WC, Towarnicki SG, Melvin RG, Youngson NA, Garvin MR, Hu Y, Nielsen S, Thomas T, Pickford R, Bustamante S, Vila-Sanjurjo A, Smyth GK, Ballard JWO (2018) Genotype to phenotype: Diet-by-mitochondrial DNA haplotype interactions drive metabolic flexibility and organismal fitness. PLoS Genet 14(11): e1007735. doi: [10.1371/journal.pgen.1007735](https://doi.org/10.1371/journal.pgen.1007735)

PROTOCOL STATUS

Working

- 1 Female third instar wandering larvae were weighed, and methanol was added to make up each sample to 20 mg/ ml.
- 2 Metabolites were extracted using an ultrasonic probe (30 s), 1 h incubation at 4Å° C and then centrifugation to remove particulates
- 3 100 Åµl aliquots of the supernatant were derivatised [200] before mass spectrometric interrogation with an Agilent GC/MSD system (Agilent Technologies, CA, USA) controlled by Chemstation software.
- 4 The GC inlet temperature was set to 230Å° C. 1.0 Åµl of derivatised sample was injected in splitless mode, using helium as a carrier gas at constant-flow of 1.0 ml/ min.
- 5 Chromatographic separation was performed on a 30 m SH-RXi-5Sil MS column (Shimadzu, NSW, Australia) with 0.25 mm internal diameter and 0.25 Åµm film thickness.
- 6 The oven temperature was programmed at 70Å° C for 2 min, then ramped at 15Å° C/ min to 320Å° C, and held 8 min.
- 7 Electron ionisation mass spectra were recorded at 1.4 scans/ s over the range m/z 50â€“700. The MSD auxiliary temperature, source temperature, and quadrupole temperature were set to 280Å° C, 230Å° C, and 150Å° C, respectively.
- 8 Analytes were identified using the NIST 2011-Wiley Mass Spectra Library. GC peaks containing mass spectra with a match quality (spectral purity) of more than 70% were considered to be tentatively identified.
- 9 No internal standard was employed. Normalisation was performed by maintaining a constant sample mass per volume.

10 Peak areas were log-transformed and statistical analysis used the limma package. Benjamini-Hochberg's correction was used to control the FDR.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited