



Mar 27,
2019

Working

Sub-Cloning Primer Design and PCR [↗](#)

Abby Jurgensmeier¹, Moriah Beck¹

¹Wichita State University

[dx.doi.org/10.17504/protocols.io.zjcf4iw](https://doi.org/10.17504/protocols.io.zjcf4iw)

Beck Lab

Abby Jurgensmeier

EXTERNAL LINK

https://www.embl.de/pepcore/pepcore_services/cloning/pcr_strategy/primer_design/amplification/index.html

PROTOCOL STATUS

Working

We use this protocol in our group and it is working

MATERIALS

NAME

CATALOG

VENDOR

dNTPs

Taq DNA Polymerase with Standard Taq Buffer - 20,000 units

M0273E

New England Biolabs

BEFORE STARTING

{ "blocks": [{ "key": "ef1qd", "text": "Primer Design for subcloning a gene of interest into the MCS of a new vector", "type": "unstyled", "depth": 0, "inlineStyleRanges": [{ "offset": 0, "length": 76, "style": "bold" }], "entityRanges": [], "data": {} }, { "key": "9pm8g", "text": "Forward", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "8n6m7", "text": "1. The forward primer overlaps with the 5'-end of the gene of interest.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "2cq7", "text": "2. Use a restriction site compatible with an available restriction enzyme near the beginning of the multiple cloning site in your vector. Add 2-10 nucleotides on the 5' end of the site for optimal cleavage efficiency.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "7ipop", "text": "3. For a C-terminal tag, include a start codon. The restriction site for the NdeI enzyme already has a start codon (CATATG), so it may be beneficial to use this site.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "1bb8e", "text": "4. Overlap with the 5' end of the gene of interest to achieve a primer of 18-30 bases total and a melting temperature (Tm) of 60°C (calculation below). The primer should have a GC content of 40-60%.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "fpnj7", "text": "", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "fpnd4", "text": "Reverse", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "dsno5", "text": "1. The reverse primer overlaps with the DNA strand complementary to the 3'-end of the gene of interest.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "e8von", "text": "2. Use a restriction site compatible with an available restriction enzyme after (3' to) the restriction site chosen above. Add 2-10 nucleotides on the 5' end of the site for optimal cleavage efficiency.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "3epma", "text": "3. For an N-terminal tag, include a stop codon. To increase termination efficiency, you may want to use 2-3 stop codons.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "809ua", "text": "4. Overlap with the 3' end of the gene of interest to achieve a primer of 18-30 bases total and a melting temperature (Tm) of 60°C (calculation below). The primer should have a GC content of 40-60%.", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "d1g0m", "text": "", "type": "unstyled", "depth": 0, "inlineStyleRanges": [], "entityRanges": [], "data": {} }, { "key": "5or2e", "text": "Primer Melting Temperature:" }] }

- 1 Combine the following in a small PCR tube:

1 uL

Forward Primer (10 uM)

1 uL	Reverse Primer (10 uM)
1uL	Template DNA (may need adjustment for at least 20 ng; remember to adjust water volume!)
1 uL	dNTPs
0.25 uL	Taq Polymerase
5 uL	10X Thermopol Buffer
41.75 uL	ddH2O

2 PCR Conditions as programmed into thermocycler as BECKTAQT

Denaturation ⌚ 00:00:30 🌡 95 °C

x30

Denaturation ⌚ 00:00:15 🌡 95 °C

Annealing ⌚ 00:00:45 🌡 55 °C **Adjust to 5 less than primer Tm**

Elongation ⌚ 00:00:30 🌡 68 °C

Final extension ⌚ 00:05:00 🌡 68 °C



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited