

RNA Extraction with Trizol for cells Version 2

Maria Cármen Sales, Maryana Branquinho, Maysa Silva

Abstract

Citation: Maria Cármen Sales, Maryana Branquinho, Maysa Silva RNA Extraction with Trizol for cells. protocols.io

dx.doi.org/10.17504/protocols.io.ssheeb6

Published: 22 Aug 2018

Before start

Clean the work area with 70% alcohol. Use all filters and autoclaved tips. Refrigerate your centrifuge to 4C.

Materials

- RNAse-free Water by Contributed by users
- TRIzol Reagent 15596026 by Thermo Fisher Scientific
- ✓ Isopropanol by Contributed by users
- Chloroform 319988 by Sigma
- ✓ 0.5% Sodium dodecyl sulfate solution by Contributed by users
- Ethanol 75% by Contributed by users
- ✓ 0.1 mM EDTA Ethylenediaminetetraacetic acid by Contributed by users

Protocol

Step 1.

Remove the culture medium from the plate and add 750 μ l of Trizol to 1 x 10⁵-10⁷ cells.

Homogenize with the pipette and transfer to a tube.

In this step, you can freeze this sample for up to 6 months in the -80C freezer or continue the protocol.

■ AMOUNT
750 μl: Trizol

I TEMPERATURE
-80 °C: freezer

Step 2.

Incubate for 5 minutes to permit complete dissociation of the nucleoproteins complex.

O DURATION

00:05:00 : Incubation

Step 3.

Add 200 µl L of chloroform per 1 mL of Trizol and incubate for 2 to 3 minutes.

O DURATION

00:02:00 : Incubation

Step 4.

Centrifuge the sample for 15 minutes at 12,000 xg at 4°C.

The mixture separates into a lower red phenol-chloroform, and interphase, and a colorless upper aqueous phase.

■ TEMPERATURE

4 °C : Centrifugation

O DURATION

00:15:00 : Centrifugation

Step 5.

Transfer the aqueous (incolor) phase containing the RNA to a new tube by angling the tube at 45° and pipetting the solution out.

Step 6.

Add 500 µl of isopropanol per 1 mL of Trizol. Incubate for 10 minutes.

O DURATION

00:10:00 : Incubation

Step 7.

Centrifuge for 10 minutes at 12,000 \times g at 4°C.

↓ TEMPERATURE

4 °C : Centrifugation

© DURATION

00:10:00 : Centrifugation

Step 8.

Discard the supernatant with a micropipetto and resuspend the pellet in 1 mL of 75% ethanol per 1 mL of Trizol.

Step 9.

Vortex the sample briefly and centrifuge for 5 minutes at 7500 \times g at 4°C.

▮ TEMPERATURE

4 °C : Centrifugation

O DURATION

00:05:00 : Centrifugation

Step 10.

Discard the supernatant with a micropipettor. Vacuum or air dry the RNA pellet for 5-10 minutes.

O DURATION

00:05:00 : Vacuum or air dry

Step 11.

Resuspend the pellet in 20–50 μ L of RNase-free water, 0.1 mM EDTA, or 0.5% SDS solution by pipetting up and down.

■ AMOUNT

20 μl : RNase-free water, 0.1 mM EDTA, or 0.5% SDS solution

Step 12.

Incubate in a water bath or heat block set at 55-60°C for 10-15 minutes.

▮ TEMPERATURE

55 °C: water bath

O DURATION

00:10:00 : water bath

Step 13.

Store in freezer -80°C until use or make the cDNA reaction then.

▮ TEMPERATURE

-80 °C: storage

Warnings

The triazole reagent is toxic, gloves, lab coat, mask should be used. Manipulate the trizol in the exhaust hood.