

Extracting Bacterial DNA from filters

Natalie Solonenko

Abstract

Adapted from Tara Oceans extraction protocol (see citation).

Citation: Natalie Solonenko Extracting Bacterial DNA from filters. **protocols.io**

dx.doi.org/10.17504/protocols.io.qv8dw9w

Published: 11 Jun 2018

Protocol

Step 1.

Add 1mL lysis buffer to each filter in a 5mL snap cap tube.

Lysis buffer:

- i. 40mM EDTA
- ii. 50mM Tris
- iii. 0.75M sucrose

Step 2.

Incubate 45min at 37°C, shaking gently

Step 3.

Add SDS to 1% v/v (in this case, 100ul of 10% SDS solution was added to each sample)

Step 4.

Incubate 1hr at 55°C, shaking gently

Step 5.

Collect liquid lysate into a 2mL tube

Step 6.

Add 1mL phenol:chloroform:isoamyl alcohol and mix well

Step 7.

Centrifuge 5min at 8000g

Step 8.

Transfer aqueous phase to new 2mL tube

Step 9.

Repeat steps 6-8

Step 10.

Add 1mL chloroform and mix well

Step 11.

Centrifuge 5min at 8000g

Step 12.

Transfer aqueous phase to upper reservoir of a 4mL 100kDa Amicon concentrator

Step 13.

Centrifuge at 1000g until <200ul remains

Step 14.

Add 2mL sterile water to upper reservoir

Step 15.

Vortex for 20sec at 1500rpm

Step 16.

Centrifuge at 1000g until <200ul remains

Step 17.

Repeat steps 14-16

Step 18.

Collect sample from upper reservoir and transfer to 1.5mL tube

Step 19.

Check concentration with Qubit HS DNA assay
