

Oct 16, 2019

Transformation of Chemically Competent Cells

NUS iGEM¹¹National University of Singapore

1

Works for me[dx.doi.org/10.17504/protocols.io.8achsaw](https://doi.org/10.17504/protocols.io.8achsaw)**NUS iGEM**
National University of Singapore ⚡

- 1 Add 1 μ l of miniprep plasmid or 5 μ l of Gibson product into competent cell

Hi-Fi Gibson Assembly
by **NUS iGEM,**
National University of Singapore

PREVIEW

RUN

- 1.1 Calculate for volumes of respective fragments to assemble based on their length and concentration

- 1.2 Add calculated volume of each fragment (maximum volume: 5 μ l)





- 1.3 Add in 5 μ l of 2x Hi-Fi DNA Assembly Master Mix into the PCR tube



NEBuilder HiFi DNA Assembly Master Mix - 10 rxns
by [New England Biolabs](#)
Catalog #: [E2621S](#)

- 1.4 Vortex to mix

- 1.5 Spin down PCR tube


- 1.6 Incubate samples at 50 °C for 00:45:00


- 2 Incubate cells in  **42 °C** water bath for  **00:00:45**
- 3 Incubate cells on ice for  **00:02:00**
- 4 Add  **1 ml** of LB media into cell sample




Preparation of LB Media
by **NUS iGEM,**
National University of Singapore


PREVIEW **RUN**

- 4.1 Weigh  **25 g** of Luria Broth Base powder.









Luria Broth Base (Miller's LB Broth Base)™, powder
by Thermo Fisher
Catalog #: 12795027


- 4.2 Add the powder into  **1 L** of water.




Water refers to sterilized deionized water

- 4.3 Autoclave entire bottle of LB media.
- 5 Incubate cells at  **37 °C** for  **01:00:00**
- 6 Spin down the cells at  **6000 x g** for  **00:01:00**
- 7 Remove  **700 µl** of the supernatant

- 8 Re-suspend pellet in remaining media
- 9 Transfer  100 µl of the culture onto the agar plate containing appropriate antibiotic







Preparation of LB Agar
by NUS iGEM,
National University of Singapore

PREVIEW

RUN



- 9.1 Weigh  32 g of LB agar powder.




LB Agar, powder (Lennox L agar)

by Thermo Fisher

Catalog #: [22700025](#)

- 9.2 Add the powder into  1 L of deionized water.




Water refers to sterilized deionized water

- 9.3 Autoclave entire bottle of LB media.

- 10 Spread the cells evenly

- 11 Incubate at  37 °C overnight

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited