

# 0.5 M EDTA (0.5 L)

Haven Himmighoefer

## Abstract

Is used in a multitude of experiments, but is often used as an ingredient in 10x TBE buffer

**Citation:** Haven Himmighoefer 0.5 M EDTA (0.5 L). **protocols.io**

dx.doi.org/10.17504/protocols.io.ks7cwhn

**Published:** 14 Nov 2017

## Before start

Make sure you are able to use the pH machine and an autoclave machine

## Materials

Ethylenediamine Tetraacetic Acid Disodium Salt Dihydrate [S311-500](#) by [Fisher Scientific](#)

Sodium Hydroxide [BP359-500](#) by [Fisher Scientific](#)

## Protocol

### Step 1.

Fill 0.5 L bottle with 300mL of deionized water.

### Step 2.

Add 95.05 g of Ethylenediamine Tetraacetic Acid Disodium Salt Dihydrate to the water and mix. The pH will be about 3 at this points. The pH needs to be at 8.

### Step 3.

Add NaOH pellets until the pH is at 8. This will be about 105 pellets. Can check if at 8 with pH machine or with pH paper.

### Step 4.

Once at pH 8, fill to 500 mL.

## **Step 5.**

Autoclave under setting 12 with water in the container to prevent over boiling.

## **Warnings**

Will be using NaOH pellets, make sure to wear gloves.