

# Quantification of circulating microRNA using TaqMan Low Density Array (TLDA)

Helle Glud Binderup, Jonna Skov Madsen, Kim Houllind, Rikke Fredslund Andersen, Claus Lohman Brasen

## Abstract

A protocol for quantification of circulating microRNA using TaqMan Low Density Array. MicroRNAs are purified from platelet-poor plasma using Nucleospin columns. During the purification step samples are spiked with cel-miR-39 as a mean of technical normalization. Normalization is performed using the  $\Delta C_t$ -method.

**Citation:** Helle Glud Binderup, Jonna Skov Madsen, Kim Houllind, Rikke Fredslund Andersen, Claus Lohman Brasen. Quantification of circulating microRNA using TaqMan Low Density Array (TLDA). **protocols.io**. dx.doi.org/10.17504/protocols.io.q62dzge

**Published:** 05 Jul 2018

## Guidelines

Blood samples should be obtained using a minimum of venous stasis and with discard of the first 3 mL of blood.

Platelet-poor plasma should be prepared within 2 hours from blood sampling.

## Before start

We recommend to test if the volume of spike-in added during microRNA purification is adequate before purifying all your samples.

When testing the preamplification product, the  $C_t$ -values of the spike-in and the target miRNA should be within the same range, otherwise adjust the volume of spike-in added during miRNA purification.

## Materials

K2-EDTA containing tubes 366643 by Becton-Dickinson

Nucleospin® miRNA Plasma 740971.50 by Macherey-nagel

TaqMan® MicroRNA Reverse Transcription Kit 4366597 by Applied Biosystems, Foster City, California, USA

2X TaqMan Universal PCR Master Mix  
4318157 by Applied Biosystems, Foster City,  
California, USA

20X TaqMan MicroRNA Assay 4440887 by  
Thermo Fisher Scientific

Custom TaqMan® Array MicroRNA Cards  
4449135 by Applied Biosystems, Foster City,  
California, USA

## Protocol

### Preparation of platelet-poor plasma (PPP)

#### Step 1.

Dual centrifugation

1. use 10 ml of EDTA anticoagulated whole blood
2. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to new tube, leaving approximately 1 mL of plasma on top of the buffy coat
4. centrifugation at 3000 g for 15 minutes (acceleration 5, brake 6, temperature 18 °C)
5. transfer plasma phase to cryo tubes, leaving approximately 1 ml of plasma in the bottom of the tube
6. store at -80 °C

#### 📌 NOTES

**Helle Glud Binderup** 22 Jun 2018

Alternatively, a prolonged single centrifugation may be used

1. use 5 ml of EDTA anticoagulated whole blood
2. centrifuge at 3000 g for 30 minutes (acceleration 5, brake 6, temperature 18 °C)
3. transfer plasma phase to cryo-tubes, leaving approximately 0.5 ml of plasma on top of the buffy coat

#### 📌 EQUIPMENT

Equipment brand:  
Hettich centrifuge  
SKU:  
4706-01  
Specifications:  
Rotina 420R

### MicroRNA purification















#### Step 2.

#### 📌 REAGENTS

Nucleospin® miRNA Plasma

**Step 3.**

Follow the instructions given by the manufacture (se notes before you start):

<b>NucleoSpin® miRNA Plasma</b>		
<b>1 Prepare sample</b>		300 µL plasma or serum* 90 µL MLP Vortex 5 s RT, 3 min
<b>2 Precipitate protein</b>	 	30 µL MPP Vortex 5 s RT, 1 min 11,000 x g, 3 min
<b>3 Transfer supernatant</b>		Transfer clear supernatant to Collection Tube (2 mL, lid)
<b>4 Adjust binding conditions</b>	 	400 µL isopropanol Vortex 5 s
<b>5 Bind RNA and DNA</b>	 	Load sample on NucleoSpin® miRNA Column RT, 2 min 11,000 x g, 30 s
<b>6 Optional: Digest DNA</b>	 	<b>Optional:</b> <b>1<sup>st</sup></b> 700 µL MW2 11,000 x g, 30 s <b>2<sup>nd</sup></b> 250 µL MW2 11,000 x g, 2 min 50 µL rDNase in Reaction Buffer for rDNase RT, 15 min
<b>7 Wash and dry</b>	 	<b>1<sup>st</sup></b> 100 µL MW1 11,000 x g, 30 s <b>2<sup>nd</sup></b> 700 µL MW2 11,000 x g, 30 s <b>3<sup>rd</sup></b> 250 µL MW2 11,000 x g, 2 min
<b>8 Elute RNA</b>	 	30 µL RNase-free H <sub>2</sub> O RT, 1 min 11,000 x g, 1 min

**⊕ NOTES**

Helle Glud Binderup 22 Jun 2018

STEP 1: use 300 µL platelet-poor plasma

STEP 4: after addition of isopropanol, add also 5 µL spike-in (cel-miR-39,  $2.75 \times 10^{-12}$  M)

STEP 6: perform the optional DNA digest

## Reverse transcription

### Step 4.



#### REAGENTS

TaqMan® MicroRNA Reverse Transcription Kit 4366597 by Applied Biosystems, Foster City, California, USA

Custom TaqMan® Array MicroRNA Cards 4449135 by Applied Biosystems, Foster City, California, USA

## Reverse transcription

### Step 5.

Prepare RT master mix:

Component	Master mix volume per 10 µL reaction*
Customs RT primer pool (10X)	1 µL
100 mM dNTPs (with dTTP)	0.27 µL
Multiscribe RT enzyme (50 U/µL)	2 µL
10x RT buffer	1 µL
MgCl <sub>2</sub> (25 mM)	1 µL
RNase inhibitor (20 U/µL)	0.1 µL
Nuclease free water	1.63 µL

\*add 10-20% excess volume

Mix gently and place on ice



#### NOTES

Helle Glud Binderup 25 Jun 2018

Remember to include a RT-negative sample (no template)  
It is also a good idea to include a RT-positive sample (a microRNA-sample included in all runs)

### Step 6.

For each RT reaction, combine 7  $\mu\text{L}$  RT master mix with 3  $\mu\text{L}$  of purified microRNA in a 0.2  $\mu\text{L}$  polypropylene reaction tube

Keep on ice

### Step 7.

Mix gently and incubate on ice for 5 minutes

### Step 8.

Transfer the reaction tubes to a thermocycler

Incubate in 40 cycles of 16  $^{\circ}\text{C}$  for 2 min, at 42  $^{\circ}\text{C}$  for 1 min and 50  $^{\circ}\text{C}$  for 1 sec.

Finish with 85  $^{\circ}\text{C}$  for 5 min and cool to 4  $^{\circ}\text{C}$ .

### Step 9.

Continue immediately to the PCR amplification or store the RT-reaction in -20  $^{\circ}\text{C}$ .

## Preamplification of cDNA

### Step 10.



#### REAGENTS

Custom TaqMan® Array  
MicroRNA Cards 4449135 by Applied  
Biosystems, Foster City, California,  
USA

TaqMan™ PreAmp Master Mix  
4391128 by Applied Biosystems,  
Foster City, California, USA

### Step 11.

Prepare preamplification reaction mix:

Component	Volume per 25 $\mu\text{L}$ reaction*
2x TaqMan PreAmp master mix	12.5 $\mu\text{L}$
Customs Primers (10X)	2.5 $\mu\text{L}$

Nuclease free water	5 $\mu$ L
<b>Total volume</b>	<b>20 <math>\mu</math>L</b>

\*add 10-20% excess volume

Mix gently

### Step 12.

Combine 20  $\mu$ L of qPCR reaction mix with 5  $\mu$ L of RT-reaction (cDNA) in a 96 well plate  
Seal plate

### Step 13.

Incubate in 95 °C for 10 min, 55 °C for 2 min and 72 °C for 2 min  
Proceed with 14 cycles of 95 °C for 15 sec and 60 °C for 4 min.  
Finish with 99.9 °C for 10 min and cool to 4 °C

Continue immediately or store at 4 °C for up to 12 hours or at -20 °C for up to 1 week.

Test preamplification product

### Step 14.

Use a specific TaqMan assay for the spike-in (Cel-miR-39, assay 000200) and one of the target miRNAs included in the Array to test the preamplification product before loading it on the array.



#### REAGENTS

20X TaqMan MicroRNA Assay  
4440887 by Thermo Fisher Scientific

2X TaqMan Universal PCR  
Master Mix 4318157 by Applied  
Biosystems, Foster City, California,  
USA

Test preamplification product

### Step 15.

Prepare qPCR reaction mix (run each sample in doublets):

Component	Volume per 20.3 $\mu$ L reaction*
-----------	--------------------------------------

2x TaqMan Universal PCR master mix	10 µL
20x TaqMan microRNA assay	1 µL
Nuclease free water	8 µL
<b>Total volume</b>	<b>19 µL</b>

\*add 10-20% excess volume

Mix gently

Test preamplification product

#### Step 16.

Combine 19 µL of qPCR reaction mix with 1.3 µL of RT-reaction (cDNA) in a 96 well plate  
Seal plate

Test preamplification product

#### Step 17.

 SOFTWARE PACKAGE

#### Sequence Detection System (SDS), 2.4

<https://www.thermofisher.com/search/results?query=sequence+analysis+software&persona=Catalog&navId=10949&refinementAction=true&focusArea=Search+All>

Test preamplification product

#### Step 18.

Transfer plate to ABI Prism 7900HT  
Use the SDS software to set up the run  
Incubate in 50 °C for 2 min and 95 °C for 10 min  
Proceed with 40 cycles of 95 °C for 15 sec and 60 °C for 60 sec.

Test preamplification product

#### Step 19.

Evaluate the Ct-values obtained

Loading of the array

#### Step 20.

Custom TaqMan® Array  
MicroRNA Cards 4449135 by Applied  
Biosystems, Foster City, California,  
USA

### Loading of the array

#### Step 21.

Prepare samples:

Component	Volume per 120 µL reaction
2x TaqMan Universal PCR II mix No AmpErase UNG	60 µL
PreAmplification product	1,2 µL
Nuclease free water	58,8 µL
<b>Total volume</b>	<b>120 µL</b>

### Loading of the array

#### Step 22.

Load 100 µL sample-dilution to each port on the array

Centrifuge The array at 1200 rpm for 2 min

Seal the array

### Real-time PCR

#### Step 23.

##### SOFTWARE PACKAGE

#### Sequence Detection System (SDS), 2.4

<https://www.thermofisher.com/search/results?query=sequence+analysis+software&persona=Catalog&navId=10949&refinementAction=true&focusarea=Search+All>

### Real-time PCR

#### Step 24.

Transfer array to ABI Prism 7900HT

Use the SDS software to set up the run

Incubate in 50 °C for 2 min and 94.5 °C for 10 min

Proceed with 40 cycles of 97 °C for 30 sec and 59.7 °C for 60 sec.



**Step 25.**

Normalization is performed using the  $\Delta\text{Ct}$ -method ( $2^{-\Delta\text{Ct}}$ )

$\Delta\text{Ct} = \text{mean Ct}_{\text{target miRNA}} - \text{mean Ct}_{\text{cel-miR-39}}$  (mean of triplets)

---