

# Long staining procedure of nuclei in *Euplotes crassus* using DAPI Version 3

Rachele Cesaroni

## Abstract

**Citation:** Rachele Cesaroni Long staining procedure of nuclei in *Euplotes crassus* using DAPI. **protocols.io**

dx.doi.org/10.17504/protocols.io.hi8b4hw

**Published:** 01 Apr 2017

## Protocol

### Step 1.

Pellet *Euplotes crassus* cells at 400 rcf for 3 minutes, and remove as much supernatant as possible by pipetting.

### ⊕ NOTES

**Rachele Cesaroni** 01 Apr 2017

Both algae and bacteria are autofluorescent. Better to have a completely starved *Euplotes crassus* culture.

### Step 2.

Add 1 ml of 2% PFA in 1X PHEM or 4% PFA in 1X PBS to the cells, and incubate them for 10 minutes at room temperature.

### Step 3.

Pellet *Euplotes crassus* cells by centrifugation at 400 rcf for 3 minutes, and remove as much supernatant as possible by pipetting.

### Step 4.

Wash cells twice with 1X PBS (400 rcf for 3 minutes each time).

### Step 5.

Add 1 ml of TBSTEM - 3% BSA and 0.5 µl of DAPI (0.1 mg/ml) to the cells, and stain for 10 minutes at room temperature.

### Step 6.

Pellet *Euplotes crassus* cells by centrifugation at 400 rcf for 3 minutes.

**Step 7.**

Add 50 µl of Prolong medium.

**Step 8.**

Place an approx. 10 µl droplet of *Euplotes crassus* cells on a slide for observation by fluorescence microscope.