

# Property Measuring spectral reflectance and transmittance (350-2500 nm) of small and/or narrow leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere

Version 3

## **Etienne Laliberté, Raymond Soffer**

#### **Abstract**

Here we describe the standardised protocol used by the Canadian Airborne Biodiversity

Observatory (CABO) to measure leaf spectral reflectance and transmittance, using an integrating sphere fitted to a portable full-range field spectroradiometer, for the special case where an individual leaf is too small and/or too narrow to entirely cover the reflectance or transmission port of the integrating sphere. Briefly, six arrays of mature, healthy and sunlit leaves from a canopy plant are arranged on a custom sample mount, and are then used for measurements of adaxial reflectance and transmittance. Reflectance measurements are referenced to a calibrated Spectralon® disk and corrected for stray light. Our leaf spectroscopy protocol is an adaptation of that of Noda et al. (2013) to the SVC DC-R/T sphere, and also builds from protocols by the Carnegie Airborne Observatory and the SVC integrating sphere user manual.

**Citation:** Etienne Laliberté,Raymond Soffer Measuring spectral reflectance and transmittance (350-2500 nm) of small and/or narrow leaves using the Spectra Vista Corporation (SVC) DC-R/T Integrating Sphere. **protocols.io** dx.doi.org/10.17504/protocols.io.q56dy9e

Published: 20 Jun 2018

#### **Guidelines**

## **Handling Spectralon®**

- **Do not touch Spectralon**® (e.g. sphere interior, reference disks, plugs) with your fingers.
- **Do not use canned air** to remove dust on the Spectralon® disk; canned air contains chemicals that can alter Spectralon®'s optical properties.
- Do not attempt to clean Spectralon® in the field, other than blowing surface dust only on the Spectralon® reference disk or sphere plugs using the Canless Air Duster System; cleaning Spectralon® requires a special procedure that should only done in the lab.

## **Equipment**

- Spectra Vista Corporation HR-1024i full-range (350-2500 nm) field spectroradiometer
- Spectra Vista Corporation 3-inch Spectralon® DC-R/T Sphere
- Semi-rugged laptop or PDA running the SVC Scan software
- Canless Air Duster System O<sub>2</sub> Hurricane (never use canned air) to remove dust from the surface of the Spectralon® reference disk
- Plastic containers with lids to temporarily store leaf arrays during measurements (optional)

## Consumables

- Nitrile gloves for handling leaves
- Whatman No. 2 filter paper (110 mm diameter)
- Acetate sheets (to make thin plastic sample mounts)
- Manila file folders (to make thin cardbboard sample platforms)
- Scotch™ Magic Tape

#### **Before start**

- 1. Consult the user manual of the spectroradiometer and the integrating sphere to set up the instrument.
- 2. The instrument should be set up in the shade, sheltered as much as possible from the elements.
- 3. All canopy plants selected for measurements should have already been tagged, identified, and georeferenced before spectroscopy measurements start.
- 4. The spectroscopist should be positionned as close as possible to the sampled plants to minimise time from collection to measurement.
- 5. The spectroscopist should be in a confortable position and have enough room around the instrument to spread leaf samples around without the risk of mixing up individual leaves during handling.
- 6. Mature, fully-developed, healthy-looking leaves from the sunlit (>3 h per day of direct sunlight) portion of the canopy are selected for spectral measurements from the bulk leaf sample (often one of a few branches). Leaves should be collected from the uppermost surface of the branch (i.e. receiving the most direct sunlight).

#### **Materials**

Manila File Folders 116723 by Staples

Write-On Transparency Film 954144 by Staples

Whatman™ Qualitative Filter Paper: Grade 2 Circles (110 mm diameter) 09-810E by Fisher Scientific

Scotch™ Magic Tape 14172 by Staples

## **Protocol**

#### Instrument Set-Up

## Step 1.

Install the integrating sphere onto the spectroradiometer.



Etienne Laliberté 04 Jun 2018

Follow the SVC integrating sphere manual p. 6-16.

#### Instrument Set-Up

## Step 2.

Position the lamp over the sphere **primary light entrance port**.

## **A** SAFETY INFORMATION

The lamp can get very hot. Grab it by the slotted heat shield.

#### **P** NOTES

Etienne Laliberté 24 Apr 2018

Make sure lamp is secured in locked position.

#### Instrument Set-Up

## Step 3.

Power the spectroradiometer on and warm up for >15 min.

#### Instrument Set-Up

#### Step 4.

Power the integrating sphere lamp and warm up for >5 min.

#### Instrument Set-Up

## Step 5.

Screw the tethered light trap on the **transmission port** sample holder.

#### NOTES

Etienne Laliberté 24 Apr 2018

The light trap can stay on the sample holder for the entire measurement session.

#### Instrument Set-Up

## Step 6.

Check lamp alignment.

## **₽** NOTES

Etienne Laliberté 04 Jun 2018

Use a thin piece of paper at the exit of the reflectance sample port (empty port) the to ensure the light beam under-fills and is centered in the reflectance port. **If it is not, then proceed to lamp alignment** as described in the SVC integrating sphere user manual, p. 29-32.

#### Software Set-Up

## Step 7.

Ensure that the spectroradiometer is connected to the computer via Bluetooth.

## **₽** NOTES

#### Etienne Laliberté 04 Jun 2018

The default Bluetooth password for SVC HR-1024i spectroradiometers is hr1024i.

#### Software Set-Up

## Step 8.

Launch SVC Scan.



## Software Set-Up

## Step 9.

Select **auto-save mode** using *File > Data Options > Auto Save.* 

#### Software Set-Up

## **Step 10.**

Select *File* > *New* to create a new working folder in your shared '**spectra**' Google Drive CABO folder.

## **P** NOTES

## Etienne Laliberté 04 Jun 2018


If you do not yet have a shared Google Drive folder for your project, contact the CABO data manager to create one (etienne.laliberte@umontreal.ca or jeremy.goimard@umontreal.ca).

Your shared CABO Google Drive folder should be named YYYYstarted-YourLastName-Degree-University, e.g. 2018-Girard-MSc-UdeM.

## Software Set-Up

#### **Step 11.**

Create a new folder named *YYYY-MM-DD-SiteID-SpectroradiometerSerialNumber* (without spaces) within that '**spectra**' folder, and keep the default file base name (= *grMMDDYY*).



#### NOTES


#### Etienne Laliberté 04 Jun 2018

The base file name is the prefix for all files that will be saved within that folder. In auto-save mode, each file will be added an incremental step, e.g., *gr053018 0000.sig*, *gr053018 0001.sig*, etc.

#### Software Set-Up

## **Step 12.**

Go to *Control* > *Setup Instrument...* to connect the spectroradiometer.



## **₽** NOTES


## Etienne Laliberté 04 Jun 2018

Different instruments may be configured under different virtual COM ports.

#### Software Set-Up

## **Step 13.**

- 1. Choose the RAW DN optic
- 2. Select 5 sec scan time
- 3. Select Auto Integration
- 4. Ensure that the date and time are correct
- 5. Check Integration Scales Raw DN Data




#### Software Set-Up

## **Step 14.**

In Control > Overlap / Matching Settings:

- 1. Select Preserve Overlapped Detector Data
- 2. Select Matching Type > None



## **Review Protocol Summary Diagram**


## **Step 15.**

Review the **document** summarising the different sphere configurations (A–I), and the scans that need to be recorded in each configuration.

#### Prepare Sample Mounts and Holders

## **Step 16.**

Build at least of six sample mounts from thin plastic (e.g. acetate sheets).



## **₽** NOTES

## Etienne Laliberté 25 Apr 2018

The sample mount is made of a thin plastic sheet with a square window that is larger than the edge of the sample port lip (e.g.  $2.5 \text{ cm} \times 2.5 \text{ cm}$ ).


## Etienne Laliberté 25 Apr 2018

Fig. 3a from Noda et al. (2013). https://doi.org/10.1111/pce.12100

#### Prepare Sample Mounts and Holders

## **Step 17.**

Build at least six few sample platforms from thin cardboard (e.g. file folders).



#### **P** NOTES

## Etienne Laliberté 25 Apr 2018

The sample platform is made of cardboard and will eventually attach to the integrating sphere to create a flat surface.


#### Etienne Laliberté 25 Apr 2018

Fig. 3c from Noda et al. (2013). https://doi.org/10.1111/pce.12100

## Prepare Leaf Arrays

#### **Step 18.**

Fix leaves onto the sample mount with tape to make leaf array #1.



#### **P** NOTES

#### Etienne Laliberté 25 Apr 2018

Fig. 3b from Noda et al. (2013). https://doi.org/10.1111/pce.12100

#### Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows.

## Prepare Leaf Arrays

## **Step 19.**

Fix leaves onto the sample mount with tape to make leaf array #2.

#### NOTES

#### Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows.

## Prepare Leaf Arrays

#### Step 20.

Fix leaves onto the sample mount with tape to make leaf array #3.

#### **P** NOTES

#### Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows.

## Prepare Leaf Arrays

#### Step 21.

Fix leaves onto the sample mount with tape to make leaf array #4.

## **₽** NOTES

#### Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows.

## Prepare Leaf Arrays

## Step 22.

Fix leaves onto the sample mount with tape to make leaf array #5.

#### **P** NOTES

## Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows.

## Prepare Leaf Arrays

## Step 23.

Fix leaves onto the sample mount with tape to make leaf array #6.

#### **P** NOTES

#### Etienne Laliberté 20 Jun 2018

Try to position as much leaf material as possible in the center of the port where the intensity of the light beam will be highest.

Leaves should be separated by about 1 mm to avoid multiple scattering among them.

If leaves are shorter than the sample mount hole (i.e. <1.25 cm long), fix leaves by their petiole with tape in one (or two) rows

## **Prepare Leaf Arrays**

## Step 24.

Fix sample mount holding leaf array #1 onto the sample platform with tape.

## **₽** NOTES

#### Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Prepare Leaf Arrays

## Step 25.

Fix sample mount holding leaf array #2 onto the sample platform with tape.

#### **P** NOTES

#### Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Prepare Leaf Arrays

## Step 26.

Fix sample mount holding leaf array #3 onto the sample platform with tape.

## **₽** NOTES

#### Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Prepare Leaf Arrays

## **Step 27.**

Fix sample mount holding leaf array #4 onto the sample platform with tape.

**₽** NOTES

Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Prepare Leaf Arrays

#### **Step 28.**

Fix sample mount holding leaf array #5 onto the sample platform with tape.

**₽** NOTES

Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Prepare Leaf Arrays

## Step 29.

Fix sample mount holding leaf array #6 onto the sample platform with tape.

**P** NOTES

Etienne Laliberté 25 Apr 2018

The sample mount holding the leaf array should be sandwiched between the two sides of the cardboard sample platform.

## Configuration A: Reflectance Mode, Leaf Array, Reference

#### Step 30.

Screw the tethered light trap on the **reflectance port** sample holder.

**₽** NOTES

Etienne Laliberté 24 Apr 2018

The light trap can stay on the sample holder for the entire measurement session.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

## **Step 31.**

Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

**P** NOTES

Etienne Laliberté 24 Apr 2018

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (= facing inside of the sphere).

#### Configuration A: Reflectance Mode Leaf Array Reference

#### **Step 32.**

Position leaf array #1 over the **transmission port** so that the **adaxial** (upper) surfaces of the leaves face into the sphere.


## Configuration A: Reflectance Mode, Leaf Array, Reference

## **Step 33.**

Collect a 'Reference Scan' in this configuration.

## Reflectance mode

Leaf array Reference radiance  $(R_{ref,a})$ 



## **P** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **reference radiance** in reflectance mode for the **leaf arrays** 1-3 ( $R_{ref,a}$ ). The reference data will be automatically saved in all successive target scan files until a new 'Reference Scan' is made.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

## **Step 34.**

Collect a 'Target Scan' in the same configuration and save the file.

## Configuration B: Reflectance Mode, Leaf Array + Paper, Reference

## **Step 35.**

Place a filter paper directly behind the sample platform holding leaf array #1.


## **₽** NOTES

Etienne Laliberté 25 Apr 2018

The leaf array #1 should remain in the same position.

## Configuration B: Reflectance Mode, Leaf Array + Paper, Reference

## **Step 36.**



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **reference radiance** in reflectance mode for the **leaf array + paper** ( $R_{ref,ap}$ ) for leaf arrays 1-3.

#### Configuration C: Reflectance Mode, Paper, Reference

## **Step 37.**

Carefully remove the sample platform holding leaf array #1 from the transmission port.


## **P** NOTES

Etienne Laliberté 25 Apr 2018

The filter paper should remain in the same position.

#### Configuration C: Reflectance Mode, Paper, Reference

## Step 38.



## **P** NOTES

Etienne Laliberté 25 Apr 2018

This corresponds to the **reference radiance** for the **filter paper** in reflectance mode ( $R_{ref,p}$ ).

## Configuration D: Reflectance Mode, Stray light

## Step 39.

Remove the filter paper from the transmission port sample holder.

#### Configuration D: Reflectance Mode, Stray light

## Step 40.

Remove the tethered calibrated Spectralon® reflectance standard from the reflectance port.


#### Configuration D: Reflectance Mode, Stray light

## **Step 41.**

Place the tethered calibrated Spectralon® reflectance standard over the **transmission port**.

## Configuration D: Reflectance Mode, Stray light

## **Step 42.**



## **₽** NOTES

Etienne Laliberté 25 Apr 2018

This corresponds to the **stray light radiance** in reflectance mode ( $R_{str}$ ).

## Configuration E: Reflectance Mode, Paper, Target

## Step 43.

Place the filter paper over the **reflectance port**.


## **₱** NOTES

Etienne Laliberté 25 Apr 2018

Focus the measurements on the same area of the paper used in previous measurements.

#### Configuration F: Reflectance Mode, Paper, Targe

## **Step 44.**



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of the **filter paper** in reflectance mode ( $R_{tar,p}$ ).


#### Configuration F: Reflectance Mode, Leaf Array #1 + Paper, Target

## Step 45.

Position the sample platform containing leaf array #1 in front the filter paper over the **reflectance port** with the **adaxial** (upper) surfaces of the leaves face into the sphere.

#### Configuration F: Reflectance Mode, Leaf Array #1 + Paper, Target

## **Step 46.**



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #1 + paper** in reflectance mode ( $R_{\text{tar,ap,1}}$ ).

## Configuration G: Reflectance Mode, Leaf Array #1, Target

## Step 47.

Remove the filter paper from the reflectance port.


#### **P** NOTES

Etienne Laliberté 04 Jun 2018

The sample platform holding leaf array #1 should remain in the same position.

## Configuration G: Reflectance Mode, Leaf Array #1, Target

## Step 48.



## **₽** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #1** in reflectance mode ( $R_{tar.a.l}$ ).

## Configuration F: Reflectance Mode, Leaf Array #2 + Paper, Target

## Step 49.

Carefully replace leaf array #1 by leaf array #2.

#### Configuration F: Reflectance Mode, Leaf Array #2 + Paper, Target


## Step 50.

Place the filter paper directly behind leaf array #2 over the **reflectance port**.

#### Configuration F: Reflectance Mode, Leaf Array #2 + Paper, Target

## **Step 51.**

Leaf array + paper Target radiance  $(R_{tar,ap,i})$ 



## **₽** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #2 + paper** in reflectance mode ( $R_{tar,ap,2}$ ).

## Configuration G: Reflectance Mode, Leaf Array #2, Target

## Step 52.

Remove the filter paper from the reflectance port.


#### NOTES

Etienne Laliberté 04 Jun 2018

The sample platform holding leaf array #2 should remain in the same position.

#### Configuration G: Reflectance Mode, Leaf Array #2, Target

## **Step 53.**



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #2** in reflectance mode ( $R_{tar,a,2}$ ).

#### Configuration F: Reflectance Mode, Leaf Array #3 + Paper, Target

## Step 54.

Carefully replace leaf array #2 by leaf array #3.

#### Configuration F: Reflectance Mode, Leaf Array #3 + Paper, Target


## Step 55.

Place the filter paper directly behind leaf array #3 over the **reflectance port**.

#### Configuration F: Reflectance Mode, Leaf Array #3 + Paper, Target

## Step 56.

Leaf array + paper Target radiance  $(R_{tar,ap,i})$ 



## **₽** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #3 + paper** in reflectance mode ( $R_{tar,ap,3}$ ).

#### Configuration G: Reflectance Mode, Leaf Array #3, Target

## Step 57.

Remove the filter paper from the reflectance port.


## **P** NOTES

Etienne Laliberté 04 Jun 2018

The sample platform holding leaf array #3 should remain in the same position.

#### Configuration G: Reflectance Mode, Leaf Array #3, Target

## **Step 58.**



## **₽** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** of **leaf array #3** in reflectance mode ( $R_{tar.a.3}$ ).

## Configuration A: Reflectance Mode, Leaf Array, Reference

#### Step 59.

Remove leaf array #3 from the reflectance port.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

## Step 60.

Remove the tethered calibrated Spectralon® reflectance standard from the **transmission port**.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

#### Step 61.

Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

#### **P** NOTES

Etienne Laliberté 24 Apr 2018

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (= facing inside of the sphere).

#### Configuration A: Reflectance Mode, Leaf Array, Reference

## Step 62.


Position leaf array #1 over the **transmission port** so that the **adaxial** (upper) surfaces of the leaves face into the sphere.

## Step 63.

Collect a 'Target Scan' in this configuration and save the file.

## Reflectance mode

Leaf array Reference radiance  $(R_{ref,a})$ 



## **P** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **reference radiance** for **leaf arrays** 4-6 in reflectance mode  $(R_{refap})$ .

## Step 64.

Place a filter paper directly behind the sample platform holding leaf array #1.


## **₽** NOTES

Etienne Laliberté 20 Jun 2018

The leaf array #1 should remain in the same position.

## **Step 65.**

Leaf array + paper Reference radiance  $(R_{\text{ref.ap}})$ 



## **P** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **reference radiance** for the **leaf array + paper** in reflectance mode ( $R_{ref,ab}$ ) for leaf arrays 4-6.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

#### Step 66.

Remove the tethered calibrated Spectralon® reflectance standard from the reflectance port.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

#### **Step 67.**

Remove leaf array #1 and filter paper from the transmission port.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

## Step 68.

Place the tethered calibrated Spectralon® reflectance standard over the **transmission port**.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

## **Step 69.**

Place the filter paper over the **reflectance port**.

## **P** NOTES

Etienne Laliberté 20 Jun 2018

Focus the measurements on the same area of the paper used in previous measurements.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

## Step 70.


Position the sample platform containing leaf array #4 in front the filter paper over the **reflectance port** with the **adaxial** (upper) surfaces of the leaves face into the sphere.

## Configuration F: Reflectance Mode, Leaf Array #4 + Paper, Target

## **Step 71.**

Collect a 'Target Scan' in this configuration and save the file.

# Reflectance mode Leaf array + paper Target radiance $(R_{tar,ap,i})$



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #4 + paper** in reflectance mode ( $R_{tar,an,4}$ ).

#### Configuration G: Reflectance Mode, Leaf Array #4, Target

## Step 72.

Remove the filter paper from the reflectance port.

## **P** NOTES

Etienne Laliberté 20 Jun 2018

The sample platform holding leaf array #4 should remain in the same position.


#### Configuration G: Reflectance Mode, Leaf Array #4, Target

## **Step 73.**

Collect a 'Target Scan' in this configuration and save the file.

#### Reflectance mode

Leaf array Target radiance  $(R_{\text{tar,a}})$ 



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #4** in reflectance mode ( $R_{tar,a,4}$ ).

## Configuration F: Reflectance Mode, Leaf Array #5 + Paper, Target

## Step 74.

Carefully replace leaf array #4 by leaf array #5.

## Configuration F: Reflectance Mode, Leaf Array #5 + Paper, Target

## Step 75.

Place the filter paper directly behind leaf array #5 over the **reflectance port**.




Etienne Laliberté 20 Jun 2018

Focus the measurements on the same area of the paper used in previous measurements.

## Configuration F: Reflectance Mode, Leaf Array #5 + Paper, Target

## Step 76.

Leaf array + paper Target radiance  $(R_{tar,ap,i})$ 



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #5 + paper** in reflectance mode ( $R_{tar,ap,5}$ ).

#### Configuration G: Reflectance Mode, Leaf Array #5, Target

## **Step 77.**

Remove the filter paper from the reflectance port.


## **₽** NOTES

Etienne Laliberté 20 Jun 2018

The sample platform holding leaf array #5 should remain in the same position.

#### Configuration G: Reflectance Mode, Leaf Array #5, Target

## Step 78.



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #5** in reflectance mode ( $R_{tar.a.5}$ ).

## Configuration F: Reflectance Mode, Leaf Array #6 + Paper, Target

## Step 79.

Carefully replace leaf array #5 by leaf array #6.

#### Configuration F: Reflectance Mode, Leaf Array #6 + Paper, Target

## **Step 80.**

Place the filter paper directly behind leaf array #6 over the **reflectance port**.

## **₽** NOTES


Etienne Laliberté 20 Jun 2018

Focus the measurements on the same area of the paper used in previous measurements.

## Configuration F: Reflectance Mode, Leaf Array #6 + Paper, Target

## Step 81.

Leaf array + paper Target radiance  $(R_{tar,ap,i})$ 



## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #6 + paper** in reflectance mode ( $R_{tar,ap,6}$ ).

#### Configuration G: Reflectance Mode, Leaf Array #6, Target

## Step 82.

Remove the filter paper from the reflectance port.

## **₱** NOTES


Etienne Laliberté 20 Jun 2018

The sample platform holding leaf array #6 should remain in the same position.

#### Configuration G: Reflectance Mode, Leaf Array #6, Target

## Step 83.

Leaf array Target radiance  $(R_{\text{tar.a}})$ 



## **P** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** of **leaf array #6** in reflectance mode ( $R_{tar.a.6}$ ).

## Configuration A: Reflectance Mode, Leaf Array, Reference

#### Step 84.

Remove leaf array #6 from the reflectance port.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

## Step 85.

Remove the tethered calibrated Spectralon® reflectance standard from the **transmission port**.

#### Configuration A: Reflectance Mode, Leaf Array, Reference

#### Step 86.

Place the tethered calibrated Spectralon® reflectance standard over the **reflectance port**.

#### **P** NOTES

## Etienne Laliberté 20 Jun 2018

Place the standard over the reflectance port so that the light beam shines directly on its reflective surface (= facing inside of the sphere).

## Configuration A: Reflectance Mode, Leaf Array, Reference

## Step 87.

Position leaf array #1 over the **transmission port** so that the **adaxial** (upper) surfaces of the leaves face into the sphere.


## Configuration A: Reflectance Mode, Leaf Array, Reference

## Step 88.

Collect a 'Target Scan' in this configuration and save the file.

## Reflectance mode

Leaf array Reference radiance ( $R_{\rm ref.a}$ )



## **P** NOTES

Etienne Laliberté 20 Jun 2018

This third reference measurement is only used to assess system stability in reflectance mode for the leaf array reference.

## Configuration B: Reflectance Mode, Leaf Array + Paper, Reference

## Step 89.

Place a filter paper directly behind the sample platform holding leaf array #1.

## **₽** NOTES


Etienne Laliberté 20 Jun 2018

The leaf array #1 should remain in the same position.

## Configuration B: Reflectance Mode, Leaf Array + Paper, Reference

## Step 90.

Leaf array + paper Reference radiance  $(R_{\text{ref.ap}})$ 



## **P** NOTES

#### Etienne Laliberté 20 Jun 2018

This third reference measurement is only used to assess system stability in reflectance mode for the leaf array + paper reference.

#### Configuration H: Transmittance Mode, Reference

## **Step 91.**

Remove leaf array #1 and filter paper from the **transmission port**.

#### Configuration H: Transmittance Mode, Reference

## Step 92.

Position the sample platform containing leaf array #1 over the **reflectance port** with the **abaxial** (lower) surface of the leaves facing into the sphere.

#### Configuration H: Transmittance Mode, Reference

## **Step 93.**

Remove the tethered calibrated Spectralon® reflectance standard from the sphere **reflectance port**.

#### Configuration H: Transmittance Mode, Reference

## Step 94.

Remove the light trap from the transmission port sample holder.

#### Configuration H: Transmittance Mode, Reference

Step 95.

Position the lamp over the sphere **transmission port**.

## **⚠** SAFETY INFORMATION

The lamp can get very hot. Grab it by the slotted heat shield.

## **₽** NOTES

Etienne Laliberté 24 Apr 2018

Make sure lamp is secured in locked position.

#### Configuration H: Transmittance Mode, Reference

## Step 96.

Install the Spectralon® plug over the **primary light port**.

## **₽** NOTES


Etienne Laliberté 24 Apr 2018

Ensure that the curved plug is placed the correct way to match the curvature of the sphere.

#### Configuration H: Transmittance Mode, Reference

## Step 97.

Collect a '**Reference Scan**' in this configuration.



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **reference radiance** in transmittance mode for the **leaf array** ( $T_{ref}$ ). The reference data will be automatically saved in all successive target scan files until a new 'Reference Scan' is made.

#### Configuration H: Transmittance Mode, Reference

## Step 98.

Collect a 'Target Scan' in this configuration and save the file.

## Configuration I: Transmittance Mode, Target

## Step 99.

Carefully remove leaf array #1 from the **reflectance port**.

## Configuration I: Transmittance Mode, Target

## Step 100.

Gently pull lamp and transmission port sample holder away from the sphere.

**A** SAFETY INFORMATION

# The lamp can get very hot. Grab it by the slotted heat shield.

## Configuration I: Transmittance Mode, Target

## Step 101.

Place the sample platform holding leaf array #1 over the **transmission port** with the abaxial (lower) surfaces of the leaves facing into the sphere.

#### Configuration I: Transmittance Mode, Target


## Step 102.

Release the transmission sample holder and move lamp back to its locked position.

## Configuration I: Transmittance Mode, Target

## Step 103.

Collect a 'Target Scan' under this configuration and save the file.



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #1 ( $T_{tar,1}$ ).

## Configuration I: Transmittance Mode, Target

## Step 104.

Carefully replace leaf array #1 by leaf array #2.

## Configuration I: Transmittance Mode, Target

## Step 105.

Collect a 'Target Scan' under this configuration and save the file.

#### NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #2 ( $T_{tar,2}$ ).

## Configuration I: Transmittance Mode, Target

## Step 106.

Carefully replace leaf array #2 by leaf array #3.

## Configuration I: Transmittance Mode, Target

## Step 107.

Collect a 'Target Scan' for leaf array #3 in this configuration and save the file.

## **₽** NOTES

Etienne Laliberté 04 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #3 ( $T_{tar,3}$ ).

## Configuration I: Transmittance Mode, Target

## Step 108.

Carefully replace leaf array #3 by leaf array #4.

## Configuration I: Transmittance Mode, Target

## Step 109.

Collect a 'Target Scan' for leaf array #4 in this configuration and save the file.

## **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #4 (T<sub>tar.4</sub>).

#### Configuration I: Transmittance Mode, Target

## Step 110.

Carefully replace leaf array #4 by leaf array #5.

## Configuration I: Transmittance Mode, Target

## Step 111.

Collect a 'Target Scan' for leaf array #5 in this configuration and save the file.

#### **P** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #5 ( $T_{tar,5}$ ).

## Configuration I: Transmittance Mode, Target

#### Step 112.

Carefully replace leaf array #5 by leaf array #6.

## Configuration I: Transmittance Mode, Target

## Step 113.

Collect a 'Target Scan' for leaf array #6 in this configuration and save the file.

#### **₽** NOTES

Etienne Laliberté 20 Jun 2018

This corresponds to the **target radiance** in **transmittance** mode for leaf array #6 ( $T_{tar,6}$ ).

#### Configuration H: Transmittance Mode, Reference

## Step 114.

Remove leaf array #6 from the transmission sample port holder.


#### Configuration H: Transmittance Mode, Reference

## Step 115.

Position the sample platform containing leaf array #1 over the **reflectance port** with the **abaxial** (lower) surface of the leaves facing into the sphere.

#### Configuration H: Transmittance Mode, Reference

## Step 116.



## **P** NOTES

Etienne Laliberté 04 Jun 2018

This second reference radiance is only used to assess system stability in transmittance mode.

#### Calculating Absolute Adaxial Reflectance of Leaf Array

## Step 117.

The equation (Noda et al. 2013; eqn. 9) for **adaxial reflectance** of leaf array i,  $\rho_{a,i}$  is

$$\rho_{a,i} = [(R_{tar,a,i} - R_{str}) \div (R_{ref,a} - R_{str})] \times \rho_{ref} \times [1 \div (1 - G_{r,i})]$$

where

 $R_{\text{tar,a},i}$  is the target radiance of leaf array i (adaxial side) in reflectance mode,

 $R_{\text{ref.a}}$  is the reference radiance used for all leaf arrays in reflectance mode,

 $R_{\rm str}$  is the stray light radiance in reflectance mode,

 $\rho_{\rm ref}$  is the absolute reflectance of the calibrated Spectralon® reflectance standard, and

 $G_{r,i}$  is the gap fraction in reflectance mode for leaf array i, which is calculated **at 400 nm** (Noda et al.

2013; eqn. 13) by

$$G_{r,i} = [((R_{tar,ap,i} - R_{str}) \div (R_{ref,ap} - R_{str})) - ((R_{tar,a,i} - R_{str}) \div (R_{ref,a} - R_{str}))] \times (\rho_{ref} \div \rho_{p})$$

where

 $R_{\text{tar,ap},i}$  is the target radiance of leaf array i + filter paper in reflectance mode,  $R_{\text{ref,ap}}$  is the reference radiance used for all leaf arrays + filter paper in reflectance mode, and  $\rho_{p}$  is the absolute reflectance of the filter paper, which is calculated (Noda et al. 2013; eqn. 3) by

$$\rho_{\rm p} = [(R_{\rm tar.p} - R_{\rm str}) \div (R_{\rm ref.p} - R_{\rm str})] \times \rho_{\rm ref}$$

where

 $R_{\text{tar,p}}$  is the target radiance of the filter paper in reflectance mode, and  $R_{\text{ref,p}}$  is the reference radiance of the filter paper in reflectance mode.

## **₽** NOTES

Etienne Laliberté 20 Jun 2018

Leaf arrays 1-3 use the first  $R_{\text{ref,a}}$  and  $R_{\text{ref,ao}}$  measurements, whereas leaf arrays 4-6 use the second measurements.

## Calculating Absolute Adaxial Transmittance of Leaf Array

## Step 118.

The equation for **adaxial transmittance** of leaf array i,  $\tau_{a,i}$  is

(modified from eqn. 15 of Noda et al. 2013 for the SVC sphere)

$$\tau_{a,i} = [ (T_{tar,i} \div T_{ref}) - G_{t,i}] \times [1 \div (1 - G_{t,i})]$$

where

 $T_{tar,i}$  is the target radiance of leaf array i in transmittance mode,

 $T_{\rm ref}$  is the reference radiance used for all leaf arrays in transmittance mode,

 $G_{t,i}$  is the gap fraction in transmittance mode for leaf array i, which is calculated **at 400 nm** (since  $\tau_{a,i} = 0$  at that wavelength) by

$$G_{t,i} = T_{tar,i} \div T_{ref}$$

## **Warnings**

The lamp of the integrating sphere can get **very hot** and should handled from its slotted base to avoid burns.