



Sep 23, 2019

## Chloroform Phenol Phage genome isolation

[Marijn Ceelen](#)<sup>1</sup><sup>1</sup>Wageningen University

1

Works for me

[dx.doi.org/10.17504/protocols.io.7kvhkw6](https://doi.org/10.17504/protocols.io.7kvhkw6)

iGEM Wageningen 2019

[Marijn Ceelen](#)  
Wageningen University

### ABSTRACT

Protocol of the isolation of genomic DNA from phage lambda and phage T7.

- 1 Use cleaned phage stock from the [phage stock](#) protocol. Add **50 µl** DNase I 10x buffer, **1 µl** DNase I (1 U/µL), and **1 µl** RNase A (10 mg/mL) for **01:30:00** at **37 °C** without shaking to remove E. coli DNA and RNA.
- 2 Add **20 µl** of **0.5 Molarity (M)** EDTA (final concentration **20 Molarity (m)**) and incubate for **00:15:00** at **75 °C** to deactivate DNase I and RNase A
- 3 Add **2 µl** Proteinase K (20 mg/mL) and **50 µl** SDS 10% to digest the phage capsid and incubate overnight at **56 °C** without shaking.
- 4 Add equal volume phenol/chloroform/iosamyl alcohol (25:24:1) and mix well (do not vortex, as genome is easily damaged)
- 5 Centrifuge at room temperature for **00:10:00** at 10000 x g
- 6 Carefully take aqueous phase
- 7 Add equal volume chloroform/isoamyl alcohol (24:1) and mix well (do not vortex, as genome is easily damaged)
- 8 Centrifuge at room temperature for **00:10:00** at 10000 x g
- 9 Carefully take aqueous phase

10 Add 1/10 volume of Sodium acetate **[M]3 Molarity (M)**

11 Add 2,5 volume of 100% Ethanol

12 Incubate at **🧊 -20 °C** for **🕒 00:30:00**


13 Centrifuge for **🕒 00:10:00** at 4668 x G

14 Discard supernatant

15 Rinse pellet with 70 % ethanol

16 Dry pellet. Heating to **🧊 50 °C** increases drying.

17 Dissolve pellet in MQ water

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited