

Fluorescent Labeling of Reovirus with Succinimidyl Ester Dyes Version 2

Bernardo Mainou

Abstract

Fluorescent labeling of reovirus virions with succinimidyl ester-conjugated fluorescent probes

Citation: Bernardo Mainou Fluorescent Labeling of Reovirus with Succinimidyl Ester Dyes. **protocols.io**

dx.doi.org/10.17504/protocols.io.kqhcv6

Published: 08 Nov 2017

Before start

You will need purified reovirus and dimethylformamide.

Sodium bicarbonate must be made on the day of labeling.

SE-conjugated dyes are not limited to Alexa dyes.

Protocol

Step 1.

Dilute 3×10^{12} reovirus particles in up to 495 μ l of fresh 0.05 M sodium bicarbonate (pH 8.5)

- Sodium Bicarbonate (FW = 84.01)
- Make a 10X Stock (0.5 M) in ddH₂O (0.42 g into 10 ml)
- Make a 1X working solution (0.05 M) in ddH₂O.

Filter Sterilize with a 0.2 μ m filter

Step 2.

Make a 1 mM stock of the succinimidyl ester Alexa fluor in dimethylformamide.

- For Alexa 546 (MW = 1159.6) (e.g. 862 μ l dimethylformamide into 1 mg of dye)

Step 3.

Add 5 μ l of 1 mM Alexa dye to virus-sodium bicarbonate solution

a. In effect, you are labeling the virus with 10 μ M of dye

Step 4.

Rock for 90 minutes in the **dark** at room temperature or 4C (we have not observed significant differences between labeling at either temperature, but 4C may enhance stability of virus)

Step 5.

Dialyze (6-8 kDa molecular weight cut-off) in PBS overnight at 4°C, replacing with the PBS after the first hour of dialysis

Step 6.

Store virus at 4°C

Note: for a Mock control, simply follow the above steps without including viral particles. In other words, add 5 μ L of 1 mM dye to 500 μ L of 0.05 M sodium bicarbonate and proceed with the labeling protocol.