



Sep 17, 2019

seqFISH Tissue Preservation V.1

Long Cai<sup>1</sup>, Nina Dar<sup>1</sup>, Yiing Lin<sup>2</sup>, Shin Lin<sup>3</sup>

<sup>1</sup>California Institute of Technology, <sup>2</sup>Washington University, Saint Louis, <sup>3</sup>University of Washington

*In Development* dx.doi.org/10.17504/protocols.io.7e5hjg6

Human BioMolecular Atlas Program (HuBMAP) Method Development Community  
Tech. support email: Jeff.spraggins@vanderbilt.edu

ABSTRACT

Tissue Preparation for seqFISH

MATERIALS TEXT


Fixation

- 1
- Clean surface with RNAzap. Make the following solutions.
- 1L PBS 1x, RNase-free
1. take 1 L RNase-free water bottle

2. use serologic pipet and draw up 100 mL water into another container

3. add 100 mL of 10x RNase-free PBS

4. mix thoroughly
- 50% sucrose
1. tare 1 L RNase-free disposable container

2. take 500 g RNase-free sucrose container

3. add RNase-free water to sucrose and transfer to 1 L container

4. continue to add RNase-free water until mass to 1 kg

5. mix thoroughly, may need agitation on platform at RT for some time

6. store at 4°C
- NOTE: for fixation, minimum PFA to tissue volume is 20:1, ideal closer to 50:1

- 2 Prepare 4% PFA in 1x PBS, 120 mL total for 4 tissues 30 mL each in 50 mL conicals.

prepare 3x conicals for 4 tissues	
16% PFA	10 ml
10x PBS	4 ml
RNAse free water	26 ml
total	40 ml

- 3 Cut freshly dissected tissue to small pieces for heart, try cutting slab, 2 cm x 1 cm x 0.4 cm thickness: 2 cm dimension from epicardium to endocardium.
- 4 Place each piece in a 50 ml tube (or smaller if you can cut to smaller pieces) in 4% PFA leave at RT for 3 hrs.  
\*\* the sample should be completely covered.
- 5 Wash three times with equal volume (30mL) of 1x PBS to remove PFA leave a bit of extra liquid in the bottom between each wash.

#### Sucrose gradient

- 6 Place the tissue in 10% Sucrose at RT and wait for it to sink in 50 mL conical, add 10 mL of sucrose wait for sample to sink to bottom – minimum incubation time of 30 min remove supernatant for next bath.
- 7 Place the tissue in 20% Sucrose at RT and wait for it to sink.
- 8 Place the tissue in 30% sucrose in 4°C – ready for shipping. Make sure there is enough 30% sucrose solution in container e.g. if in 50 mL conicals, have 50 mL of 30% sucrose.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited